

Dimensions:  
Width: 8-3/4" (22,2 cm)  
Length: 11-1/8" (28,3 cm)

**PB-10001-000A**  
**PBZ-10001-000A**

T-12703-100C  
Width: 8.7598 in = 222.5mm  
Length: 10.374 in = 263.5mm

**ARROW**



**REF PW-18080**

(IPNXXXXXX)

**.018** inch dia.  
spring-wire  
guide


**80** cm  
spring-wire  
guide

**Spring-Wire Guide .018" (.46 mm) Diameter x 31-1/2" (80 cm) Wire Proximal End: Solid Wire .018" Diameter Distal End: Coil .018" Diameter x 7 cm**

**Contents:**

1: Spring-Wire Guide: .018" (0.46 mm) dia. x 31-1/2" (80 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other)

**Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU**  
Fluid path components are non-pyrogenic.



EU Authorized  
Representative and  
Importer:

**EC REP** Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
2797 Co. Westmeath, Ireland

Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**LOT** Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(11)YYMMDD  
(10)Lot Number

Product of xxx  
Packaged in xxx

**Teleflex**  
www.teleflex.com

LBL077152 R00 (2025-09)



**REF PW-18080**

**.018** inch dia.  
spring-wire  
guide

**80** cm  
spring-wire  
guide

**ARROW**

**Spring-Wire Guide .018" (.46 mm) Diameter x 31-1/2" (80 cm) Wire Proximal End: Solid Wire .018" Diameter Distal End: Coil .018" Diameter x 7 cm**