

Dimensions: 6-3/8" (16 cm) x 8-1/4" (21 cm)



REF RC-09700 (IPNXXXXXX)

7 Fr. 5.08 cm catheter length .025 inch dia. spring-wire guide

RIC® Rapid Infusion Catheter Exchange Set



Contents:

- 1: Sheath: 7 Fr. x 2" (5.08 cm) Arrow-Flex® Radiopaque Polyurethane with Tissue Dilator
1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)
1: Scalpel: #11

All components are CE 2797 unless otherwise noted.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
Not made with natural rubber latex. Fluid path components are non-pyrogenic.



EU Authorized Representative and Importer:

Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
2797 Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Product of xxx
Packaged in xxx

Teleflex®
www.teleflex.com

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

LBL077154 R00 (2025-08)