

Dimensions:
 Width: 6-9/64" (15.6 cm)
 Length: 8-21/32" (22 cm)

PB-50014-100B
PBZ-50014-103A

T50014100A6 Dimensions: 5 7/8" x 7 5/8"



REF CV-16702-E
 (IPNXXXXXX)



Two-Lumen CVC

Contents:


- 1: Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXX
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX
- 1: Pressure Transduction Probe Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 2: Dust Cap: Non-Vented Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: SecondSite™ Adjustable Hub: Fastener Origin: XX XXXX XX XXX XXXXXXXX XXXXXXXXXXXX
- 1: SecondSite™ Adjustable Hub: Catheter Clamp Origin: XX XXXX XX XXX XXXXXXXX
XXXXXXXXXXXX

All components are CE 2797 unless otherwise noted.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU California Prop. 65

⚠ WARNING: Cancer and Reproductive Harm.
 www.P65Warnings.ca.gov

Not made with natural rubber latex. Fluid path components are non-pyrogenic.



Priming
 Volume*
 (mL)

Gravity
 Flow
 Rate†
 (mL/hr)

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr)
Distal (16 Ga.)	0.47	2681
Proximal (16 Ga.)	0.52	3013

* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.

CAS # 7440-48-4

EU Authorized Representative and Importer:

EC REP Teleflex Medical
 IDA Business and
 Technology Park
 Dublin Road, Athlone
 2797 Co. Westmeath, Ireland

All components are CE 2797 unless otherwise noted.

Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA

Set Product of xxx
 Package in xxx

Teleflex®
 www.teleflex.com

LOT Lot Number

YYYY-MM-DD

YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXXX
 (17)YYMMDD
 (11)YYMMDD
 (10)Lot Number

LBL080150 R00 (2026-02)



7 Fr. 2 Lumen

20 cm catheter length

.032 inch dia. spring-wire guide

REF CV-16702-E

Two-Lumen CVC

