



REF **SAC-01618-CA**

(IPNXXXXXX)

18 Ga.

16 cm catheter length

.025 inch dia. spring-wire guide

Arterial Catheterization



Contents:

- 1: Catheter: 18 Ga. x 16 cm Radiopaque PEBA
- 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-3/4" (7 cm) XTW
- 1: Injection Needle: Eclipse™¹ 25 Ga. x 1" (2.54 cm)
- 1: Injection Needle: SafetyGlide™¹ 23 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: Syringe: 5 mL Luer-Slip
- 1: SharpsAway® Disposal Cup
- 1: SharpsAway® II Locking Disposal Cup
- 1: Drape: 41" x 55" with 4" fenestration
- 1: Towel: 13" x 18" (33 cm x 46 cm)
- 1: Needle Holder
- 1: Safety Scalpel: #11
- 1: Checklist Sheet
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: Tegaderm®² 10 cm x 12 cm
- 1: Surgical Apparel: Impervious Gown
- 1: Surgical Apparel: Bouffant Cap
- 1: Surgical Apparel: Mask with Eye Shield
- 1: Suture: 3-0 Silk with Curved Needle

¹A trademark of Becton, Dickinson and Company.

Catheter Priming Volume: 0.24 mL Priming volume is approximate.

Not made with natural rubber latex.

