

Dimensions: 6-3/8" (16 cm) x 8-1/4" (21 cm)

ARROW®

REF RC-09850

(IPNXXXXXX)

8.5 Fr.

6.4 cm
catheter
length.025 inch dia.
spring-wire
guide

RIC® Rapid Infusion Catheter Exchange Set

en fr it no es

Contents:

- 1: Sheath: 8.5 Fr. x 2-1/2" (6.4 cm) Arrow-Flex® Radiopaque FEP with Tissue Dilator
1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)
1: Scalpel: #11

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. www.teleflex.com/IFU
Not made with natural rubber latex. Fluid path components are non-pyrogenic.



Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Product of xxx
Packaged in xxx

Teleflex®
www.teleflex.com

LOT Lot Number



YYYY-MM-DD



YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(11)YYMMDD
(10)Lot Number

LBL080684 R00 (2025-09)