




# Arrow® Pressure Injectable Multi-Lumen PICC Kit

REF EU-05563-HPMSB

**Contents:**

- 1: Multi-Lumen Indwelling Catheter: 6 Fr. x 55 cm **Pressure Injectable** Radiopaque Polyurethane with Blue FlexTip®, T-Port Connector, Extension Line Clamps, Contamination Guard and Placement Wire
- 1: Peel-Away Sheath: 13 Ga. x 7 cm Radiopaque over 6.5 Fr. Dilator
- 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 33 cm (Straight Soft Tip on Both Ends)
- 1: Catheter: 22 Ga. x 4.45 cm Radiopaque over 25 Ga. RW Introducer Needle
- 1: Filter Needle: 18 Ga. x 3.81 cm TW with Blunt Fill Tip
- 1: **Injection Needle: Eclipse™ 22 Ga. x 3.81 cm RW**
- 1: Introducer Needle: 21 Ga. x 3.81 cm TW
- 1: Introducer Needle: Echogenic 21 Ga. x 7 cm TW
- 1: Syringe: 10 mL Luer-Lock
- 1: Syringe: 10 mL Luer-Slip
- 1: **Injection Needle: Eclipse™ 25 Ga. x 2.54 cm RW and 3 mL Luer-Lock Syringe**
- 1: Foil-Laminated Pouch 79% Skin Protectant Prep Pad
- 3: Dust Caps: Non-Vented
- 1: Clamp: Catheter
- 1: Fastener: Catheter Clamp
- 1: **SharpsAway II™ Locking Disposal Cup**
- 1: **Catheter Trimmer**
- 2: CSR Wraps
- 1: **Drape: 91 cm x 104 cm with 8 cm x 13 cm fenestration**
- 1: **Drape: 152 cm x 193 cm**
- 1: Towel
- 1: Towel: 61 cm x 91 cm
- 1: **Safety Scalpel: #11**
- 1: Chart Sticker<sup>5</sup>
- 1: Flow Rate Information Card with Injection Log<sup>5</sup>
- 1: Patient ID Card<sup>5</sup>
- 1: Patient Information Booklet<sup>5</sup>
- 1: PICC Line Insertion Checklist<sup>5</sup>
- 1: Sterile Procedure Sign<sup>5</sup>
- 1: **Dressing: STATLOCK® Catheter Stabilization Device**
- 1: Dressing: Tegaderm®<sup>4</sup> 15.5 cm x 10 cm
- 5: Gauze Pads: 5 cm x 5 cm
- 10: Gauze Pads: 10 cm x 10 cm
- 2: Paper Tape Measures (one packaged externally)
- 1: **Surgical Apparel: Bouffant Cap**
- 1: **Surgical Apparel: Impervious Gown**
- 1: **Surgical Apparel: Mask**<sup>5</sup>
- 1: **Surgical Apparel: Mask with Eye Shield**
- 1: Tourniquet<sup>5</sup>
- 1: Tape: Steri-Strip®<sup>4</sup>
- 1: HemoHopper® Fluid Receptacle

<sup>1</sup> A trademark of Becton Dickinson and Company.  
<sup>2</sup> US Patent Nos. 6,123,193 and 6,276,527.  
<sup>3</sup> A registered trademark of C. R. Bard, Inc.  
<sup>4</sup> A registered trademark of 3M Company.  
<sup>5</sup> Packaged externally to the sterile kit.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	Pump Flow Rate‡ (mL/hr)	MAX Pressure Injection Flow Rate** (mL/sec)
Distal (17 Ga.)	0.71	1070	5976	6
Medial (19 Ga.)	0.46	360	2246	not rated
Proximal (19 Ga.)	0.49	390	2467	not rated

\* Priming volumes are approximate and are done without injection caps.  
 † Flow rates were determined using room temperature water, 100 cm head height and represent approximate flow capabilities.  
 ‡ Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.  
 \*\* Pressure injection flow rates are determined at the injector pressure setting of 300 psi maximum using media of 11.8 centipoise viscosity, with 152 cm pressure tubing.



Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

- NL** Kit met PICC met meerdere lumina voor hogedrukinjectie, met Blue FlexTip katheter, positioneringsdraad, prik- en snijbeveiliging en voorzorgsmaatregelen voor maximale steriele barrière
- FR** Kit PICC multi-lumière pour injection sous pression avec cathéter Blue FlexTip, fil guide, accessoires de protection contre les objets piquants, coupants, tranchants et « précautions barrières maximales »
- DE** Mehrlumiges Pressure Injectable PICC (Peripherally Inserted Central Catheter) Kit mit Blue FlexTip Katheter, Platzierungsdraht, Nadel- und Klängen-Entsorgungsbehälter mit Verschluss und maximalen Barriere-Vorsichtsmaßnahmen
- EL** Kit PICC πολλαπλών αυλών για έγχυση υπό πίεση, που περιλαμβάνει καθετήρα με άκρο Blue FlexTip, σύρμα τοποθέτησης, διατάξεις προστασίας από αιχμηρά αντικείμενα και μέσα διαφύλαξης μέγιστων στειρών προφυλάξεων
- IT** Kit con PICC a lume multiplo per infusione a pressione con catetere Blue FlexTip, funzioni di sicurezza per taglienti e precauzioni per massima barriera di sterilità
- PL** Zestaw cewnika PICC nadającego się do wstrzykiwania pod ciśnieniem, wielokanałowego, z końcówką Blue FlexTip, drutem umieszczającym, zabezpieczeniami przed ukłuciem ostrymi narzędziami i środkami ostrożności zapewnianymi maksymalnie sterylą osłonę
- PT** Kit de cateter central de inserção periférica de lúmenes múltiplos para injeção a pressão com cateter com ponta Blue FlexTip, estilete de introdução, funcionalidades de segurança para objectos cortantes e precauções de protecção integral
- RU** Многопросветный комплект PICC для введения веществ под давлением с катетером Blue FlexTip, размещителем, предохранительным устройством для острых предметов и предупредительными мерами максимального барьера
- ES** Kit de cateterización central de inserción periférica multiluz para inyección a presión con catéter con punta Blue FlexTip, fiador para colocación, características de seguridad frente a objetos punzantes y precauciones de máxima barrera
- SV** Set med tryckinjicerbar perifert införd central kateter med flera lumen och Blue FlexTip, placeringsstråd, säkerhetsfunktioner för stickskydd och försiktighetsåtgärder för maximal barriär
- TR** Basınçla Enjekte Edilebilir Çok Lümenli PICC Kit, Blue FlexTip Kateterli, Yerleştirme Telli, Kesici Madde Güvenlik Özellikli ve Maksimum Bariyer Önlemleri



EC REP Teleflex Medical  
 IDA Business and Technology Park  
 Athlone, Ireland



P-05563-100B (12/15)



## Arrow® Pressure Injectable Multi-Lumen PICC Kit

REF EU-05563-HPMSB

**6** Fr. | **3** Lumen | **55** cm catheter length | **.018** inch dia. spring-wire guide

