



Large-Bore Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter for High Volume Infusions

REF **CS-16123-F**

12 Fr.

3 Lumen

25 cm catheter length

.035 inch dia. spring-wire guide

Contents:


- 1: Multi-Lumen Indwelling Catheter: 12 Fr. x 25 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
- 3: Dust Cap: Non-Vented
- 1: Drape: 24" x 36" (61 cm x 91 cm)
- 1: Scalpel: #11

All components are CE 0086 unless otherwise noted.

Rx only.

Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.41	2530
Medial (12 Ga.)	1.50	11190
Proximal (12 Ga.)	1.50	12400

* Priming volumes are approximate and are done without accessories.
 † Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

- fr** Jeu de cathétérisme veineux central de grand calibre à lumières multiples avec cathéter Blue FlexTip pour les perfusions de grand volume
- de** Besteck für viellumigen, zentralen, großkalibrigen Venenkatheter mit Blue-FlexTip-Katheter für Infusionen großer Flüssigkeitsmengen
- it** Set per cateterismo della vena centrale multilume di diametro ampio con catetere Blue FlexTip per volumi di infusione elevati
- pl** Wielokanałowy zestaw o dużej średnicy do cewnikowania żył centralnych z cewnikiem Blue FlexTip do przetaczania dużych objętości płynów infuzyjnych
- pt** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Grande Calibre, com Cateter Blue FlexTip, para Infusões de Grandes Volumes
- ru** Комплект для катетеризации главных вен с многопросветным катетером большого внутреннего диаметра с Blue FlexTip для вливаний большого объема
- sl** Komplet za centralno vensko katetrizacijo z velikimi luknjami in več svetlinami s katetrom Blue FlexTip za infundiranje velike prostornine
- es** Conjunto de cateterización venosa central de luz múltiple y gran diámetro con catéter Blue FlexTip para infusiones de gran volumen
- sv** Set för storhålig flerkanalig central venkatetrering med Blue FlexTip-kateter för infusioner med hög volym
- tr** Büyük Açıklıklı Multi-Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile, Yüksek Hacimli İnfüzyonlar için



Not made with natural rubber latex

Manufacturer:
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA

EC REP Teleflex Medical
 IDA Business and Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland
 0086

Teleflex®
 www.teleflex.com

S-16123-132D (07/19)



Large-Bore Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter for High Volume Infusions

REF **CS-16123-F**

12 Fr.

3 Lumen

25 cm catheter length

.035 inch dia. spring-wire guide