



## Central Venous Catheter with Blue FlexTip® for use only with Arrow MAC™ Two-Lumen Central Venous Access Device

REF **ML-00701**

**7** Fr. | **14** cm indwelling catheter length with MAC™


### Contents:

- 1: Indwelling Catheter: 7 Fr. Radiopaque Polyurethane with Blue FlexTip®, Catheter Contamination Guard, Extension Line Clamps
  - 1: Dust Cap: Non-Vented
- All components are CE 0086 unless otherwise noted.

Rx only.

**Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (12 Ga.)	0.66	9310

\* Priming volumes are approximate and are done without accessories.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

- fr** Cathéter veineux central avec cathéter Blue FlexTip à utiliser exclusivement avec le dispositif d'accès veineux central à deux lumières Arrow MAC
- de** Zentraler Venenkatheter mit Blue FlexTip nur zum Gebrauch mit dem Einführbesteck für den zweilumigen, zentralen Venenkatheter von Arrow MAC
- it** Catetere per vena centrale con Blue FlexTip da usare esclusivamente con il dispositivo di accesso venoso centrale a due lumi Arrow MAC
- pl** Cewnik do cewnikowania żył centralnych z końcówką Blue FlexTip, do użytku tylko z dwukanałowym urządzeniem Arrow MAC do uzyskiwania dostępu do żył centralnych
- pt** Cateter Venoso Central com Blue FlexTip para utilizar apenas com o Dispositivo para Acesso Venoso Central com Lúmen Duplo Arrow MAC
- ru** Катетер для главных вен с Blue FlexTip, для применения исключительно с двухпросветным устройством для доступа в главные вены Arrow MAC
- sl** Centralni venški kateter z Blue FlexTip je samo za uporabo s pripomočkom za venški dostop Arrow MAC z dvema svetlinama
- es** Catéter para vena central con Blue FlexTip para uso exclusivo con el dispositivo de acceso venoso central de dos luces de Arrow MAC
- sv** Centralvenös kateter med Blue FlexTip, endast för användning med Arrow MAC tvåkanalig produkt för centralt ventillträde
- tr** Santral Venöz Kateter, Blue FlexTip ile, sadece Arrow MAC İki Lümenli Santral Venöz Erişim Cihazı ile kullanılmak içindir

**CE** REP 0086  
Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

**Manufacturer:**  
Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex®**  
www.teleflex.com

IL-00701-130C (2020-02)



## Central Venous Catheter with Blue FlexTip® for use only with Arrow MAC™ Two-Lumen Central Venous Access Device

REF **ML-00701**

**7** Fr. | **14** cm indwelling catheter length with MAC™