

ARROW®



Insertion Set for Arrow 4.5 Fr. PICCs

REF **EU-00745-45**

Contents:

- 1: GlideThru™ Peel-Away Sheath: 4.5 Fr. x 7 cm Radiopaque over 4.5 Fr. Dilator
 - 1: Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
 - 1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
 - 1: SharpsAway II™ Locking Disposal Cup
 - 1: Safety Scalpel: #11
- nl** Inbengset voor Arrow 4,5 Fr. PICC's
 - fr** Kit d'insertion pour PICC 4,5 Fr. Arrow
 - de** Einführset für 4,5 Fr. Arrow PICCs
 - it** Set di inserimento per Arrow 4,5 Fr. PICCs
 - pl** Zestaw wprowadzający do cewników PICC 4,5 Fr. firmy Arrow
 - pt** Conjunto de Inserção para Cateteres Centrais de Inserção Periférica (PICCs) de Calibre 4,5 Fr. da Arrow
 - ru** Набор для введения PICCs Arrow 4,5 Fr.
 - es** Set de inserción para PICC de 4,5 Fr. de Arrow
 - sv** Införingsset för Arrow 4,5 Fr. PICCs
 - tr** Arrow 4,5 Fr. PICCs için Yerleştirme Seti



Not made with
natural rubber latex

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
CE
0086

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

SZ-00745-101B (07/18)



Insertion Set for Arrow 4.5 Fr. PICCs
REF **EU-00745-45**

ARROW®