

Dimensions:
 Width: 10-5/8" (27 cm)
 Length: 12-1/2" (31.8 cm)

Contents:

- 1: Midline Catheter, Trimmable: 4 Fr x 8" (20 cm)
 Radiopaque Polyurethane with T-Part Connector,
 Integral Extension Line, Extension Line Clamp,
 Integral Suture Wing, Contamination Guard
 and Placement Wire
- 1: Spring-Wire Guide, Nitinol, Marked: .018"
 (0.46 mm) dia. x 17-3/4" (45 cm) (Straight
 Soft Tip on One End - Straight Stiff Tip on Other)
 with Arrow Advancer™
- 1: Introducer Needle: Echogenic 21 Ga x 2-3/4"
 (7.00 cm) TW
- 1: Syringe: 10 ml Luer-Lock
- 1: GlideThru™ Peel-Away Sheath: 4 Fr x 2-3/4"
 (7 cm) Radiopaque over 4 Fr Dilator
- 1: Dressing: GRIP-LOK® Catheter Fixation Device
- 1: Safety Scalpel: #11
- 1: Catheter Trimmer
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Touriquet
- 1: Paper Tape Measure
- 1: ... Valve: MicroCLAVE®² Neutral Displacement Connector

¹ A registered trademark of TDI.
² A registered trademark of ICU Medical, Inc.

Fluid path components are non-pyrogenic.



Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA

EC REP Teleflex Medical
 IDA Business and Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland

Teleflex®
 www.teleflex.com

LOT




REF EU-02041-ML

Arrow® Midline Kit

ARROW®

- NL** Arrow midline katheterkit
- FR** Kit Arrow Midline
- DE** Arrow Midline-Kit
- IT** Kit catetere Midline Arrow
- PT** Kit de Linha Mediana da Arrow
- ES** Kit de catéter de línea media de Arrow
- SV** Arrow Midline-sats



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	Pump Flow Rate‡ (mL/hr)
Distal (18 Ga)	0.30	1690	7110

* Priming volumes are approximate and are done without injection caps.
 † Flow rates were determined using room temperature water, 100 cm head height and represent approximate flow capabilities.
 ‡ Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.
 Above values apply to the untrimmed catheter. Trimming catheter will alter priming volumes and flow rates.

SE02041-0108 (12/15)



REF

EU-02041-ML

Arrow® Midline Kit

4 Fr

20 cm catheter length

.018 inch dia. spring-wire guide

