

Dimensions:  
Width: 6-3/4" (17,2 cm)  
Length: 11-5/16" (28,7 cm)

REF EU-04150-N

## Central Venous Catheterization Set

# ARROW®

### Contents:

- 1: ..... Indwelling Catheter: 20 Ga x 5" (12 cm) Radiopaque Polyurethane with Integral Suture Wing
- 1: ..... Spring-Wire Guide, Nitinol, Marked: .025" (0.64 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer™
- 1: ..... Introducer Needle: 20 Ga x 1-1/2" (4 cm) XTW
- 1: ..... Syringe: 3 ml Luer-Slip
- 1: ..... Tissue Dilator: 4 Fr (1.35 mm) x 7.6 cm
- 1: ..... Tissue Dilator: 4 Fr (1.35 mm) x 4.1 cm
- 1: ..... SecondSite™ Adjustable Hub: Fastener
- 1: ..... SecondSite™ Adjustable Hub: Catheter Clamp
- 1: ..... Safety Scalpel: #11
- 1: ..... SharpsAway II™ Locking Disposal Cup

All components are CE 0086 unless otherwise noted.

Fluid path components are non-pyrogenic.



Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

EC REP Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

**Teleflex®**  
www.teleflex.com

- BE** Набор для центральной венозной катетеризации
- HR** Centralni venski kateter Komplet
- CS** Centrální žilní katetrizační set
- DA** Centralvenøs kateteriserings-sæt
- NL** Centraal veneuze katheterisatieset
- ET** Tsentraalveeni kateteriseerimiskomplekt
- FI** Keskuslaskimokatetri sarja
- FR** Jeu de cathétérisme veineux central
- DE** Besteck für zentralen Venenkatheter
- EL** Σετ Καθετηριασμού Κεντρικής φλέβας
- HU** Centrális vénás katéterezési szett
- IS** Holæðaleggsett
- IT** Set per cateterismo della vena centrale
- LT** Centrinės venos kateterizavimo rinkinys
- NO** Sentralt venekateteriseringssett
- PL** Zestaw do cewnikowania żył centralnych
- PT** Conjunto de Cateterização Venosa Central
- RO** Set de cateterizare venoasă centrală
- RU** Комплект для установки центрального венозного катетера
- SR** Komplet za Centralnu vensku kateterizaciju
- SK** Centrálna venózna katetrizačná súprava
- SL** Centralni venski komplet za katetrizacijo
- ES** Conjunto de cateterización venosa central
- SV** Set för central venkateterisering
- TR** Santral Venöz Kateterizasyon Seti
- UK** Комплект для встановлення центрального венозного катетера

Catheter Diameter: 20 Ga, 1.1 mm

Priming Volume: 0.12 ml

Priming volumes are approximate and are done without the injection cap.

Flow Rate: 1900 ml/hr

Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

**SHARPS  
SAFETY  
FEATURES**

5704150-1048 (12/15)

LOT



ARROW  
INTERNATIONAL



REF EU-04150-N

20 Ga

12 cm  
catheter  
length

.025 inch dia.  
spring-wire  
guide

Central Venous Catheterization Set

Peel to  
Open