

Dimensions:
Width: 6-3/4" (17,2 cm)
Length: 11-5/16" (28,7 cm)

REF **EU-04301-EN**

Central Venous Catheterization Set with Blue FlexTip® Catheter

ARROW®

Contents:

- 1: Indwelling Catheter: 16 Ga x 8" (20 cm) Radiopaque Polyurethane with Blue FlexTip®, Integral Extension Line, Extension Line Clamp, Integral Suture Wing
- 1: Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End — "J" Tip on Other) with Arrow Advancer™ with ECG Mark
- 1: Introducer Needle: 18 Ga x 2-1/2" (6.35 cm) XTW
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 ml
- 1: Tissue Dilator: 7 Fr (2.3 mm) x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Pressure Transduction Probe

All components are CE 0086 unless otherwise noted.

Fluid path components are non-pyrogenic.



 Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

 Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

 www.teleflex.com

- be** Набор для центральной венозной катетеризации с катетером Blue FlexTip
- hr** Centralni venski kateter Komplet s Blue FlexTip kateter
- cs** Centrální žilní katetrizační set s katetrem Blue FlexTip
- da** Centralvenøs kateteriserings-sæt med Blue FlexTip kateter
- nl** Centraal veneuze katheterisatieset met Blue FlexTip katheter
- et** Tsentraalveeni kateteriseerimiskomplekt Blue FlexTip kateetriga
- fi** Keskuslaskimokatetrisarja Blue FlexTip katetrilla
- fr** Jeu de cathétérisme veineux central avec cathéter Blue FlexTip
- de** Besteck für zentralen Venenkatheter mit Blue FlexTip Katheter
- el** Σετ Καθητηριασμού Κεντρικής φλέβας με Καθετήρα Blue FlexTip
- hu** Centrális vénás katéterezési szett, Blue FlexTip katéterrel
- is** Holæðaleggsett með Blue FlexTip æðalegg
- it** Set per cateterismo della vena centrale con catetere Blue FlexTip
- lt** Centrinės venos kateterizavimo rinkinys su „Blue FlexTip“ kateteriu
- no** Sentralt venekateteriseringssett med Blue FlexTip kateter
- pl** Zestaw do cewnikowania żył centralnych z cewnikiem z końcówką Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Cateter Blue FlexTip
- ro** Set de cateterizare venoasă centrală cu cateter Blue FlexTip
- ru** Комплект для установки центрального венозного катетера с катетером Blue FlexTip
- sr** Komplet za Centralnu vensku kateterizaciju sa Blue FlexTip kateterom
- sk** Centrálna venózna katetrizačná súprava s katétrom Blue FlexTip
- sl** Centralni venski komplet za katetrizacijo s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central con catéter Blue FlexTip
- sv** Set för central venkateterisering med Blue FlexTip kateter
- tr** Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- uk** Комплект для встановлення центрального венозного катетера з катетером Blue FlexTip

Catheter Diameter: 16 Ga, 1.7 mm

Priming Volume: 0.43 ml

Priming volumes are approximate and are done without the injection cap.

Flow Rate: 2900 ml/hr

Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

91/010 8611-1090725

LOT



ARROW



REF

EU-04301-EN

16^{Ga} 20^{cm} catheter length .032^{inch dia.} spring-wire guide

Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open