

Dimensions:
Width: 6-3/4" (17.2 cm)
Length: 11-5/16" (28.7 cm)

REF **EU-04522-N**

Central Venous Catheterization Set

ARROW®

Contents:

- 1: Indwelling Catheter: 22 Ga x 4" (10 cm) Radiopaque Polyurethane with Integral Extension Line, Extension Line Clamp, Integral Suture Wing
- 1: Spring-Wire Guide, Nitinol, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends) with Arrow Advancer™
- 1: Introducer Needle: 20 Ga x 1-1/2" (4 cm) XTW
- 1: Syringe: 5 ml Luer-Slip
- 1: Tissue Dilator 4 Fr (1.35 mm) x 1-5/8" (4.1 cm)
- 1: Tissue Dilator 4 Fr (1.35 mm) x 3" (7.6 cm)
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter
- 1: Safety Scalpel: #11
- 1: SharpsAway II™ Locking Disposal Cup

All components are CE 0086 unless otherwise noted.

Contains no medication.
Fluid path components are non-pyrogenic.



 Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

 Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

- BE** Набор для центральной венозной катетеризации
- HR** Centralni venski kateter Komplet
- CS** Centrální žilní katetrizační set
- DA** Centralvenøs kateteriserings-sæt
- NL** Centraal veneuze katheterisatieset
- ET** Tsentraalveeni kateteriseerimiskomplekt
- FI** Keskuslaskimokatetrisarja
- FR** Jeu de cathétérisme veineux central
- DE** Besteck für zentralen Venenkatheter
- EL** Σετ Καθητηριασμού Κεντρικής φλέβας
- HU** Centrális vénás katéterezési szett
- IS** Holæðaleggsett
- IT** Set per cateterismo della vena centrale
- LT** Centrinės venos kateterizavimo rinkinys
- NO** Sentralt venekateteriseringssett
- PL** Zestaw do cewnikowania żył centralnych
- PT** Conjunto de Cateterização Venosa Central
- RO** Set de cateterizare venoasă centrală
- RU** Комплект для установки центрального венозного катетера
- SR** Komplet za Centralnu vensku kateterizaciju
- SK** Centrálna venózna katetrizačná súprava
- SL** Centralni venski komplet za katetrizacijo
- ES** Conjunto de cateterización venosa central
- SV** Set för central venkateterisering
- TR** Santral Venöz Kateterizasyon Seti
- UK** Комплект для встановлення центрального венозного катетера

**SHARPS
SAFETY
FEATURES**

Catheter Diameter: 22 Ga, 0.9 mm

Priming Volume: 0.21 ml

Priming volumes are approximate and are done without the injection cap.
Injection cap priming volume is 0.17 ml.

Flow Rate: 1300 ml/hour

Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

52-04522-100A (1/14)

LOT



(01) 0 08 01902 08766 9

ARROW
INTERNATIONAL



REF
EU-04522-N

Central Venous Catheterization Set

22^{Ga}

10^{cm}
catheter
length

.021

inch dia.
spring-wire
guide

Peel to
Open