

Dimensions:  
Width: 6-3/4" (17,2 cm)  
Length: 11-5/16" (28,7 cm)

REF **EU-04701-N**

## Central Venous Catheterization Set with Blue FlexTip® Catheter

# ARROW®

### Contents:

- 1: ..... Indwelling Catheter: 14 Ga x 8" (20 cm) Radiopaque Polyurethane with Blue FlexTip®, Integral Extension Line, Extension Line Clamp, Side Holes, Integral Suture Wing
- 1: ..... Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End — "J" Tip on Other) with Arrow Advancer™ with ECG Mark
- 1: ..... Introducer Needle: 18 Ga x 2-1/2" (6.35 cm) XTW
- 1: ..... Syringe: 5 ml Luer-Slip
- 1: ..... Tissue Dilator: 7 Fr (2.3 mm) x 10.2 cm
- 1: ..... SecondSite™ Adjustable Hub: Fastener
- 1: ..... SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 0086 unless otherwise noted.

Fluid path components are non-pyrogenic.



 Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

 Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

 **Teleflex®**  
www.teleflex.com

- be** Набор для центральной венозной катетеризации с катетером Blue FlexTip
- hr** Centralni venski kateter Komplet s Blue FlexTip kateter
- cs** Centrální žilní katetizační set s katetrem Blue FlexTip
- da** Centralvenøs kateteriserings-sæt med Blue FlexTip kateter
- nl** Centraal veneuze katheterisatieset met Blue FlexTip katheter
- et** Tsentraalveeni kateteriseerimiskomplekt Blue FlexTip kateetriga
- fi** Keskuslaskimokatetrisarja Blue FlexTip katetrilla
- fr** Jeu de cathétérisme veineux central avec cathéter Blue FlexTip
- de** Besteck für zentralen Venenkatheter mit Blue FlexTip Katheter
- el** Σετ Καθετηριασμού Κεντρικής φλέβας με Καθετήρα Blue FlexTip
- hu** Centrális vénás katéterezési szett, Blue FlexTip katéterrel
- is** Holæðaleggsett með Blue FlexTip æðalegg
- it** Set per cateterismo della vena centrale con catetere Blue FlexTip
- lt** Centrinės venos kateterizavimo rinkinys su „Blue FlexTip“ kateteriu
- no** Sentralt venekateteriseringssett med Blue FlexTip kateter
- pl** Zestaw do cewnikowania żył centralnych z cewnikiem z końcówką Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Cateter Blue FlexTip
- ro** Set de cateterizare venoasă centrală cu cateter Blue FlexTip
- ru** Комплект для установки центрального венозного катетера с катетером Blue FlexTip
- sr** Komplet za Centralnu vensku kateterizaciju sa Blue FlexTip kateterom
- sk** Centrálna venózna katetizačná súprava s katétrom Blue FlexTip
- sl** Centralni venski komplet za katetrizacijo s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central con catéter Blue FlexTip
- sv** Set för central venkatetrering med Blue FlexTip kateter
- tr** Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- uk** Комплект для встановлення центрального венозного катетера з катетером Blue FlexTip

Catheter Diameter: 14 Ga, 2.2 mm

Priming Volume: 0.55 ml

Priming volumes are approximate and are done without the injection cap.

Flow Rate: 4900 ml/hr

Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

LOT



SC4701-1208 (10/16)

**ARROW**



REF

EU-04701-N

14<sup>Ga</sup>

20<sup>cm</sup>  
catheter  
length

.032<sup>inch dia.</sup>  
spring-wire  
guide

Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to  
Open