

Dimensions:  
Width: 8-3/4" (22.2 cm)  
Length: 11-1/8" (28.3 cm)

**REF EU-04701-CVTS**  
**Central Venous Catheterization Set**  
**with Blue FlexTip® Catheter**

# ARROW®

**Contents:**

- 1: .....Indwelling Catheter: 14 Ga x 8" (20 cm) Radiopaque Polyurethane with Blue FlexTip®, Integral Extension Line, Extension Line Clamp, Side Holes, Integral Suture Wing
- 1: .....Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™ with ECG Mark
- 1: .....Introducer Needle: Echogenic 18 Ga x 2-1/2" (6.35 cm) XTW
- 1: .....Syringe: 5 ml Luer-Slip
- 1: .....Tissue Dilator: 7 Fr (2.3 mm) x 10.2 cm
- 1: .....Fastener: Catheter Clamp
- 1: .....Clamp: Catheter
- 1: .....Safety Scalpel: #11
- 1: .....SharpsAway II™ Locking Disposal Cup
- 1: .....ECG Extension Cable

All components are CE 0086 unless otherwise noted.

Fluid path components are non-pyrogenic. Contains no medication.



Arrow International, Inc.  
2400 Bernville Road  
Reading, PA 19605 USA

**EC REP** Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

- BE** Набор для центральной венозной катетеризации з катетэрам Blue FlexTip
- HR** Centralni venski kateter Komplet s Blue FlexTip kateter
- CS** Centrální žilní katetrizační set s katetrem Blue FlexTip
- DA** Centralvenøs kateteriserings-sæt med Blue FlexTip kateter
- NL** Centraal veneuze katheterisatieset met Blue FlexTip katheter
- ET** Tsentraalveeni kateteriseerimiskomplekt Blue FlexTip kateetriga
- FI** Keskuslaskimokatetri sarja Blue FlexTip katetrilla
- FR** Jeu de cathétérisme veineux central avec cathéter Blue FlexTip
- DE** Besteck für zentralen Venenkatheter mit Blue FlexTip Katheter
- EL** Σετ Καθητηριασμού Κεντρικής φλέβας με Καθετήρα Blue FlexTip
- HU** Centrális vénás katéterezési szett, Blue FlexTip katéterrel
- IS** Holæðaleggsett með Blue FlexTip æðalegg
- IT** Set per cateterismo della vena centrale con catetere Blue FlexTip
- LT** Centrinės venos kateterizavimo rinkinys su „Blue FlexTip“ kateteriu
- NO** Sentralt venekateteriseringssett med Blue FlexTip kateter
- PL** Zestaw do cewnikowania żył centralnych z cewnikiem z końcówką Blue FlexTip
- PT** Conjunto de Cateterização Venosa Central com Cateter Blue FlexTip
- RO** Set de cateterizare venoasă centrală cu cateter Blue FlexTip
- RU** Комплект для установки центрального венозного катетера с катетером Blue FlexTip
- SR** Komplet za Centralnu vensku kateterizaciju sa Blue FlexTip kateterom
- SK** Centrálna venózna katetrizačná súprava s katétrom Blue FlexTip
- SL** Centralni venski komplet za katetrizacijo s katetrom Blue FlexTip
- ES** Conjunto de cateterización venosa central con catéter Blue FlexTip
- SV** Set för central venkatetrering med Blue FlexTip kateter
- TR** Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- UK** Комплект для встановлення центрального венозного катетера з катетером Blue FlexTip

Catheter Diameter: 14 Ga, 2.2 mm

Priming Volume: 0.55 ml

Priming volumes are approximate and are done without the injection cap.

Flow Rate: 4900 ml/hour

Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

**LOT**



(01) 0 08 01902 07349 5

52-04701-126A (12/14)



**REF** EU-04701-CVTS

Central Venous Catheterization Set with Blue FlexTip® Catheter

14<sup>Ga</sup> 20<sup>cm</sup> catheter length

.032<sup>inch dia.</sup> spring-wire guide

Peel to Open