

Dimensions:
 Width: 6-3/4" (17,2 cm)
 Length: 11-5/16" (28,7 cm)

PMS 193

REF ES-04706-P

**Central Venous
 Catheterization Set
 with Blue FlexTip® Catheter**

ARROW®

Contents:

- 1: Indwelling Catheter: 14 Ga x 6" (16 cm) Radiopaque Polyurethane with Blue FlexTip®, Integral Extension Line, Extension Line Clamp, Side Holes, Integral Suture Wing
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End — "J" Tip on Other) with Arrow Advancer™
- 1: Introducer Needle: 18 Ga x 2-1/2" (6.35 cm) XTW and 5 ml Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Tissue Dilator: 7 Fr (2.3 mm) x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Pressure Transduction Probe

All components are CE 0086 unless otherwise noted.

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

Sterile. Fluid path components are non-pyrogenic.

Manufacturer
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA

CE **0086**
ECIREP Teleflex Medical
 IDA Business and Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland

Teleflex®
 www.teleflex.com

- fr** Jeu de cathétérisme veineux central avec cathéter Blue FlexTip
- de** Besteck für zentralen Venenkatheter mit Blue FlexTip Katheter
- it** Set per cateterismo della vena centrale con catetere Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Cateter Blue FlexTip
- es** Conjunto de cateterización venosa central con catéter Blue FlexTip
- sv** Set för central venkatetrising med Blue FlexTip kateter

Catheter Diameter: 14 Ga, 2.2 mm

Priming Volume: 0.48 ml

Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.

Flow Rate: 5400 ml/hour

Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

PMS 193

ES-04706-103F(5/16)



Caution



Consult instructions for use



Do not reuse



Do not resterilize

STERILE EO

Sterilized by ethylene oxide



Keep away from sunlight



Keep dry



Do not use if package is damaged



Not made with natural rubber latex

LOT

Use By

PMS 193

PMS 193

ARROW®
 INTERNATIONAL



REF ES-04706-P

14^{Ga}

16 cm
 catheter length

.032 inch dia.
 spring-wire guide

Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open