

Dimensions:
 Width: 6-3/4" (17,2 cm)
 Length: 11-5/16" (28,7 cm)

REF EU-04730-CVT

**Central Venous
 Catheterization Set**
 with Blue FlexTip® Catheter

ARROW®

Contents:

- 1: Indwelling Catheter: 14 Ga x 12" (30 cm) Radiopaque Polyurethane with Blue FlexTip®, Integral Extension Line, Extension Line Clamp, Side Holes, Integral Suture Wing
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™ with ECG Mark
- 1: Introducer Needle: Echogenic 18 Ga x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 ml Luer-Slip
- 1: Tissue Dilator: 7 Fr (2.3 mm) x 10.2 cm
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter
- 1: Scalpel: #11
- 1: ECG Extension Cable

All components are CE 0086 unless otherwise noted.

Contains no medication.
 Fluid path components are non-pyrogenic.



Arrow International, Inc.
 2400 Bernville Road
 Reading, PA 19605 USA

EC REP Teleflex Medical
 IDA Business and Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland

- BE** Набор для центральной венозной катетеризации з катэтрам Blue FlexTip
- HR** Centralni venski kateter Komplet s Blue FlexTip kateter
- CS** Centrální žilní katetrizační set s katetrem Blue FlexTip
- DA** Centralvenøs kateteriserings-sæt med Blue FlexTip kateter
- NL** Centraal veneuze katheterisatieset met Blue FlexTip katheter
- ET** Tsentraalveeni kateteriseerimiskomplekt Blue FlexTip katetriga
- FI** Keskuslaskimokatetri sarja Blue FlexTip katetrilla
- FR** Jeu de cathétérisme veineux central avec cathéter Blue FlexTip
- DE** Besteck für zentralen Venenkatheter mit Blue FlexTip Katheter
- EL** Σετ Καθετηριασμού Κεντρικής φλέβας με Καθετήρα Blue FlexTip
- HU** Centrális vénás katéterezési szett, Blue FlexTip katéterrel
- IS** Holæðaleggset með Blue FlexTip æðalegg
- IT** Set per cateterismo della vena centrale con catetere Blue FlexTip
- LT** Centrinės venos kateterizavimo rinkinys su „Blue FlexTip“ kateteriu
- NO** Sentralt venekateteriseringssett med Blue FlexTip kateter
- PL** Zestaw do cewnikowania żył centralnych z cewnikiem z końcówką Blue FlexTip
- PT** Conjunto de Cateterização Venosa Central com Cateter Blue FlexTip
- RO** Set de cateterizare venoasă centrală cu cateter Blue FlexTip
- RU** Комплект для установки центрального венозного катетера с катетером Blue FlexTip
- SR** Komplet za Centralnu vensku kateterizaciju sa Blue FlexTip kateterom
- SK** Centrálna venózna katetrizačná súprava s katétrom Blue FlexTip
- SL** Centralni venski komplet za katetrizacijo s katetrom Blue FlexTip
- ES** Conjunto de cateterización venosa central con catéter Blue FlexTip
- SV** Set för central venkateterisering med Blue FlexTip kateter
- TR** Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- UK** Комплект для встановлення центрального венозного катетера з катетером Blue FlexTip

Catheter Diameter: 14 Ga, 2.2 mm

Priming Volume: 0.75 ml
 Priming volumes are approximate and are done without the injection cap.

Flow Rate: 5300 ml/hour
 Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

574047301308 (10/14)

LOT



(01) 0 08 01902 06572 8

ARROW
 INTERNATIONAL



REF

EU-04730-CVT

14 Ga

30 cm catheter length

.032 inch dia. spring-wire guide

Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open