

Dimensions:
Width: 6-3/4" (17,2 cm)
Length: 11-5/16" (28,7 cm)

REF **EU-04730-N**

Central Venous Catheterization Set with Blue FlexTip® Catheter

ARROW®

Contents:

- 1: Indwelling Catheter: 14 Ga x 12"
(30 cm) Radiopaque Polyurethane
with Blue FlexTip®, Integral Extension Line,
Extension Line Clamp, Side Holes,
Integral Suture Wing
- 1: Spring-Wire Guide, Nitinol, Marked:
.032" (0.81 mm) dia. x 26-3/4"
(68 cm) (Straight Soft Tip on One End -
"J" Tip on Other) with Arrow Advancer™
with ECG mark
- 1: Introducer Needle:
18 Ga x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 ml Luer-Slip
- 1: Tissue Dilator: 7 Fr (2.3 mm)
x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub:
Catheter Clamp

All components are CE 0086 unless
otherwise noted.

Fluid path components are non-pyrogenic.



 Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

 Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Teleflex®
www.teleflex.com

LOT



- be** Набор для центральной венозной катетеризации с катетером Blue FlexTip
- hr** Centralni venski kateter Komplet s Blue FlexTip kateter
- cs** Centrální žilní katetrizační set s katetrem Blue FlexTip
- da** Centralvenøs kateteriserings-sæt med Blue FlexTip kateter
- nl** Centraal veneuze katheterisatieset met Blue FlexTip katheter
- et** Tsentraalveeni kateteriseerimiskomplekt Blue FlexTip kateetriga
- fi** Keskuslaskimokatetrisarja Blue FlexTip katetrilla
- fr** Jeu de cathétérisme veineux central avec cathéter Blue FlexTip
- de** Besteck für zentralen Venenkatheter mit Blue FlexTip Katheter
- el** Σετ Καθετηριασμού Κεντρικής φλέβας με Καθετήρα Blue FlexTip
- hu** Centrális vénás katéterezési szett, Blue FlexTip katéterrel
- is** Holæðaleggsett með Blue FlexTip æðalegg
- it** Set per cateterismo della vena centrale con catetere Blue FlexTip
- lt** Centrinės venos kateterizavimo rinkinys su „Blue FlexTip“ kateteriu
- no** Sentralt venekateteriseringssett med Blue FlexTip kateter
- pl** Zestaw do cewnikowania żył centralnych z cewnikiem z końcówką Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Cateter Blue FlexTip
- ro** Set de cateterizare venoasă centrală cu cateter Blue FlexTip
- ru** Комплект для установки центрального венозного катетера с катетером Blue FlexTip
- sr** Komplet za Centralnu vensku kateterizaciju sa Blue FlexTip kateterom
- sk** Centrálna venózna katetrizačná súprava s katétrom Blue FlexTip
- sl** Centralni venski komplet za katetrizacijo s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central con catéter Blue FlexTip
- sv** Set för central venkatetrering med Blue FlexTip kateter
- tr** Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- uk** Комплект для встановлення центрального венозного катетера з катетером Blue FlexTip

Catheter Diameter: 14 Ga, 2.2 mm

Priming Volume: 0.75 ml

Priming volumes are approximate and are done without the injection cap.

Flow Rate: 5300 ml/hr

Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

SL047301348 (10/16)

ARROW®



REF

EU-04730-N

14 Ga

30 cm
catheter
length

.032 inch dia.
spring-wire
guide

Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to
Open