

REF SI-09875-E

Percutaneous Sheath Introducer Set

with Integral Hemostasis Valve/Side Port
for use with 7.5 Fr. Catheters

ARROW®

Contents:


- 1: Sheath: 8.5 Fr. x 4" (10 cm) Arrow-Flex®
Radiopaque Polyurethane with Integral
Hemostasis Valve/Side Port and
Tissue Dilator
- 1: Spring-Wire Guide, Marked:
.035" (0.89 mm) dia. x 17-3/4" (45 m)
(Straight Soft Tip on One End –
"J" Tip on Other) with Arrow Advancer™
- 1: Introducer Needle:
18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire
Introduction Syringe: 5 mL
- 1: Cath-Gard® Catheter Contamination Shield:
80 cm with TwistLock™ Adapter
- 1: Obturator: 8 Fr.
- 3: Gauze Pads: 4" x 4"

Rx only.


Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Single use – do not resterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains no medication.

This Product Contains Dry Natural Rubber.

 Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

 Teleflex Medical
IDA Business and Technology Park
Athlone, Ireland



- CS** Set s percutánním zaváděcím pouzdem a integrovaným hemostatickým ventilem a bočním portem pro použití s katetry 7,5 Fr.
- DA** Perkutan sheath introducer-sæt med integreret hæmostaseventil/sideåbning til brug med 7,5 Fr. katetre
- NL** Percutane hulsintroducerset met integrale hemostaseklep/zijpoort voor gebruik met 7,5 Fr. katheters
- FR** Jeu d'introduction de gaines par voie percutanée avec valve hémostatique/port latéral intégrés pour emploi avec les cathéters 7,5 Fr.
- DE** Perkutanes Schleuseneinführbesteck mit eingebauter Hämostaseventil/Seitenanschluß-Baugruppe zum Gebrauch mit Kathetern der Größen 7,5 Fr.
- EL** Σετ διαδερμικού εισαγωγέα θηκάρσιου με ενιαία αμοστατική βαλβίδα/πλευρική θύρα για χρήση με καθετήρες 7,5 Fr.
- HU** Percutan hüvelybevezető felszerelés egybeépített vérzéscsillapító szeleppel és oldalcsatlakozóval 7,5 Fr. katéterekkel való használatra
- IT** Set per introduzione di guaina percutanea con valvola emostatica/bocchetta laterale incorporate per l'uso con cateteri da 7,5 Fr.

- PL** Zestaw do przezskórnego wprowadzania koszulki, z integralnym zaworem hemostaticznym/portem bocznym do użycia z cewnikami o rozmiarze 7,5 F.
- PT** Conjunto de Introdução de Bainha Percutânea com Dispositivo de Válvula de Hemostase Integral e Porta Lateral, para utilizar com Cateteres 7,5 Fr.
- RU** Комплект чрескожной капсулы-интубатора с несъемным гемостазным клапаном/боковым портом, для использования с катетерами 7,5 Fr.
- SK** Perkutánná sada zavádzača puzdra s integrálnym hemostatickým ventilom/postranným portom pre použitie s 7,5 Fr. katétami
- SL** Uvajalni komplet s percutanim tulcem z vgrajenim hemostatskim ventilom/stranskim vhodom za uporabo s katetri 7,5 Fr.
- ES** Conjunto introductor de vaina percutánea con válvula hemostática y orificio lateral integrales para uso con catéteres de 7,5 Fr.
- SV** Set för perkutan skyddsintroducerare med integral hemostasventil/sidoöppning för bruk med katetrar mellan 7,5 Fr.
- TR** Perkütan Kılıf İntroducer Seti, Entegre Hemostaz Valfi/Yan Portlu, 7,5 Fr. Kateterler ile kullanılmak üzere

1109875-1308, Rev. 1 (6/12)

ARROW®

INTERNATIONAL

REF SI-09875-E 8.5^{Fr.} 10^{cm} sheath length 3.6^{mm} sheath O.D. .035^{inch} dia. spring-wire guide

Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port

Peel to Open