

Special Order Product



REF MTO-09903-KU

9^{Fr.} | 8.5^{cm sheath length} | 3.9^{mm sheath O.D.} | .025^{inch dia. spring-wire guide}

Percutaneous Sheath Introducer Kit
with Integral Hemostasis Valve/Side Port for use with 7 - 8 Fr. Catheters

Contents:

- 1: Sheath: 9 Fr. x 3-3/8" (8.5 cm) Arrow-Flex® Radiopaque Polyurethane with Integral Hemostasis Valve/Side Port and Tissue Dilator
- 1: Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle
- 1: Injection Needle: 25 Ga. x 1-1/2" (3.81 cm)
- 1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW
- 1: Syringe: 10 mL Luer-Slip Loss of Resistance, Plastic
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Obturator: 8 Fr.
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
- 1: Medicine Cup - clear
- 2: Drape: 24" x 36" (61 cm x 91 cm)
- 1: Needle Holder
- 1: Scalpel: #11
- 3: Stopcock: 4-way High-Flow
- 1: Dressing: Adhesive 5" x 4"(12.5 cm x 10 cm)
- 2: Prep Sponge Swab
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 2-0 Silk with Curved Needle
- 1: Tray: Prep

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.
Not made with natural rubber latex.
Fluid path components are non-pyrogenic.

FOR EXPORT ONLY



LOT



Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA
Teleflex®
www.teleflex.com
SZ-09903-128B (11/18)



REF MTO-09903-KU

9^{Fr.} | 8.5^{cm sheath length} | 3.9^{mm sheath O.D.} | .025^{inch dia. spring-wire guide}

Percutaneous Sheath Introducer Kit
with Integral Hemostasis Valve/Side Port for use with 7 - 8 Fr. Catheters

