

Dimensions:  
 Width: 6-3/4" (17,2 cm)  
 Length: 11-5/16" (28,7 cm)

PMS 259

**REF BR-10703**

**Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter**

**ARROW®**

**Contents:**

- 1: ..... Multi-Lumen Indwelling Catheter: 7 Fr x 4-3/8" (11 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: ..... Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™
- 1: ..... Introducer Needle: 18 Ga x 2-1/2" (6.35 cm) XTW and 5 ml Luer-Slip Syringe
- 1: ..... Tissue Dilator: 8.5 Fr (2.8 mm) x 10.2 cm
- 1: ..... SecondSite™ Adjustable Hub: Fastener
- 1: ..... SecondSite™ Adjustable Hub: Catheter Clamp
- 1: ..... Scalpel: #11
- 3: ..... Stopcock: 4-way High-Flow

All components are CE 0086 unless otherwise noted.

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.


Sterile. Fluid path components are non-pyrogenic.

PMS 259

**Manufacturer**  
 Arrow International, Inc.  
 Subsidiary of Teleflex Incorporated  
 2400 Bernville Road  
 Reading, PA 19605 USA

**CE REP**  
 Teleflex Medical  
 IDA Business and Technology Park  
 Dublin Road, Athlone  
 Co. Westmeath, Ireland

**Teleflex®**  
 www.teleflex.com



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga)	0.34	3600
Medial (18 Ga)	0.33	2200
Proximal (18 Ga)	0.35	2300

\*Priming volume are approximate and done without the stopcock.  
 †Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

(91/1) 9001-0201-175



Caution



Consult instructions for use



Do not reuse



Do not resterilize



Sterilized by ethylene oxide



Keep away from sunlight



Keep dry



Do not use if package is damaged



Not made with natural rubber latex

**LOT**

**Use By**

PMS 259

PMS 259

**ARROW**  
 INTERNATIONAL



**REF BR-10703**

**7 Fr 3 Lumen 11 cm catheter length .032 inch dia. spring-wire guide**

**Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter**

Peel to Open