

Dimensions:  
 Width: 6-3/4" (17,2 cm)  
 Length: 11-5/16" (28,7 cm)

PMS 259

REF **BR-10854**

**Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter**

**ARROW®**

**Contents:**

- 1: Quad-Lumen Indwelling Catheter: 8.5 Fr x 4-3/8" (11 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™
- 1: Catheter: 18 Ga x 2-1/2" (6.35 cm) Radiopaque over 20 Ga Introducer Needle, RW
- 1: Introducer Needle: 18 Ga x 2-1/2" (6.35 cm) XTW and 5 ml Luer-Slip Syringe
- 1: Tissue Dilator: 9 Fr (3.0 mm) x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Scalpel: #11
- 4: Stopcock: 4-way High-Flow

All components are CE 0086 unless otherwise noted.

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.


Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

Sterile. Fluid path components are non-pyrogenic.

**Manufacturer**  
 Arrow International, Inc.  
 Subsidiary of Teleflex Incorporated  
 2400 Bernville Road  
 Reading, PA 19605 USA

**CE** REF  
 Teleflex Medical  
 IDA Business and Technology Park  
 Dublin Road, Athlone  
 Co. Westmeath, Ireland

**Teleflex®**  
 www.teleflex.com



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga)	0.35	5300
Medial 1 (14 Ga)	0.39	8400
Medial 2 (18 Ga)	0.33	2900
Proximal (18 Ga)	0.37	2800

\*Priming volume are approximate and done without the stopcock.  
 †Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

- fr** Jeu de cathétérisme veineux central à quatre lumières avec cathéter Blue FlexTip
- de** Besteck für vierlumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- it** Set per cateterismo della vena centrale a quattro lumi con catetere Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Lúmen Quádruplo e Cateter Blue FlexTip
- es** Conjunto de cateterización venosa central de cuatro luces con catéter Blue FlexTip
- sv** Set för fyrkanalig central venkatetrisering med Blue FlexTip kateter

PMS 259

(91/5100E) 5100E 75



Caution



Consult instructions for use



Do not reuse



Do not resterilize



Sterilized by ethylene oxide



Keep away from sunlight



Keep dry



Do not use if package is damaged



Not made with natural rubber latex

LOT

Use By

PMS 259

PMS 259

**ARROW**  
 INTERNATIONAL



REF

**BR-10854**

**8.5<sup>Fr</sup> 4 Lumen**

**11 cm**  
 catheter length

**.032** inch dia.  
 spring-wire guide

Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open