

ARROW®



Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

REF **CV-12702-BZ** | 7 Fr. | 2 Lumen | 16 cm catheter length | **.025** inch dia. spring-wire guide

Contents:

- 1: Two-Lumen Indwelling Catheter: 7 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

- pt** Conjunto de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip
- es** Conjunto de cateterización venosa central de dos lúmenes con catéter Blue FlexTip


All components are CE 0086 unless otherwise noted.

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (18 Ga.)	0.34	1470
Proximal (14 Ga.)	0.51	5480

* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
0086

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

IV-12702-131B (11/19)



Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

REF **CV-12702-BZ**

7 Fr. | 2 Lumen | 16 cm catheter length | **.025** inch dia. spring-wire guide

ARROW®