

Dimensions:
Width: 6-3/4" (17.2 cm)
Length: 11-5/16" (28.7 cm)

REF EU-12703-N

Multi-Lumen Central Venous
Catheterization Set
with Blue FlexTip® Catheter

ARROW®

Contents:

- 1: Multi-Lumen Indwelling Catheter:
7 Fr x 6" (16 cm) Radiopaque
Polyurethane with Blue FlexTip®,
Extension Line Clamps
- 1: Spring-Wire Guide, Nitinol,
Marked: .032" (0.81 mm) dia.
x 17-3/4" (45 cm) (Straight Soft
Tip on One End - "J" Tip on Other)
with Arrow Advancer™
- 1: Introducer Needle: 18 Ga
x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 ml Luer-Slip
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter
- 1: Tissue Dilator 8.5 Fr (2.8 mm)
x 4" (10.2 cm)

All components are CE 0086 unless otherwise noted.


Fluid path components are non-pyrogenic.
Contains no medication.



Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

- BE** Многопросветный набор для центральной венозной катетеризации с катетером Blue FlexTip
- HR** Višelumenski centralni venski kateter Set s Blue FlexTip kateterom
- CS** Vícelumenový centrální žilní katetrizační set s katetrem Blue FlexTip
- DA** Multilumen centralvenos kateteriseringsset med Blue FlexTip-kateter
- NL** Multi-lumen centraal veneuze katheterisatieset met Blue FlexTip katheter
- ET** Mitme valendikuga tsentraalse veenikateeteri komplekt koos Blue FlexTip katetriga
- FI** Monilumeninen keskuslaskimokatetrisarja jossa Blue FlexTip-katetri
- FR** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- DE** Besteck für viellumigen, zentralen Venenkatheter mit Blue-FlexTip-Katheter
- EL** Σετ Καθετηριασμού Κεντρικής φλέβας Πολλαπλούς Αυλούς με Καθετήρα Blue FlexTip
- HU** Többlumenű, centrális vénás katéterezési szett, Blue FlexTip katéterrel
- IS** Fjöllúmena holæðaleggset með Blue FlexTip æðalegg
- IT** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- LT** Daugiakanalis centrinės venos kateterizavimo mini rinkinys su „Blue FlexTip“ kateteriu
- NO** Multilumen sentralt venekateteriseringssett med Blue FlexTip kateter
- PL** Zestaw do cewnikowania żyły centralnej z cewnikami wielokanalowym z końcówką Blue FlexTip
- PT** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip
- RO** Set de cateterizare venoasă centrală multi-lumen cu cateter Blue FlexTip
- RU** Набор с многопросветным центральным венозным катетером с кончиком Blue FlexTip
- SR** Višelumenski centralni venski kateterizacioni Set sa Blue FlexTip kateterom
- SK** Viacľúmenová centrálna venózna katetrizačná súprava s katetrom Blue FlexTip
- SL** Večlumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip
- ES** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- SV** Set för flerkanalig central venkateterisering med Blue FlexTip-kateter
- TR** Multi Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- UK** Багатопроевітний центральний венозний катетер Набор: катетер Blue с кінчиком Blue FlexTip Catheter



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)
Distal (16 Ga)	0.39	3400
Medial (18 Ga)	0.35	1800
Proximal (18 Ga)	0.37	1900

*Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.

†Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

SZ-12703-19A (12/12)

LOT



(01) 0 08 01902 07940 4



REF EU-12703-N

7 Fr 3 Lumen 16 cm catheter length .032 inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open