

Dimensions:
Width: 6-3/4" (17.2 cm)
Length: 11-5/16" (28.7 cm)

REF EU-12703-EN

Multi-Lumen Central Venous Catheterization Set

with Blue FlexTip® Catheter

ARROW®

Contents:

- 1: Multi-Lumen Indwelling Catheter: 7 Fr x 6" (16 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer™
- 1: Introducer Needle: 18 Ga x 2-1/2" (6.35 cm) XTW
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 ml
- 1: Pressure Transduction Probe
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter
- 1: Tissue Dilator 8.5 Fr (2.8 mm) x 4" (10.2 cm)

All components are CE 0086 unless otherwise noted.


Fluid path components are non-pyrogenic. Contains no medication.



Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

- CS** Vícélumenový centrální žilní katetrizační set s katetrem Blue FlexTip
- DA** Multilumen centralvenøs kateteriseringsæt med Blue FlexTip-kateter
- NL** Multi-lumen centraal veneuze katheterisatieset met Blue FlexTip katheter
- FR** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- DE** Besteck für viellumigen, zentralen Venenkatheter mit Blue-FlexTip-Katheter
- EL** Σετ Καθετηριασμού Κεντρικής φλέβας Πολλαπλούς Αυλούς με Καθετήρα Blue FlexTip
- IT** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- NO** Multilumen sentralt venekateteriseringssett med Blue FlexTip kateter
- PL** Zestaw do cewnikowania żyły centralnej z cewnikiem wielokanałowym z końcówką Blue FlexTip
- PT** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip
- RU** Набор с многопросветным центральным венозным катетером с кончиком Blue FlexTip
- SL** Večlumenski centralni venksi komplet za katetrizacijo s katetrom Blue FlexTip
- ES** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- SV** Set för flerkanalig central venkatetrering med Blue FlexTip-kateter
- TR** Multi Lúmen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)
Distal (16 Ga)	0.39	3400
Medial (18 Ga)	0.35	1800
Proximal (18 Ga)	0.37	1900

*Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.

†Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

SZ-12703-121B (6/13)

LOT



(01) 0 08 01902 07939 8



REF EU-12703-EN 7 Fr 3 Lumen 16 cm catheter length .032 inch dia. spring-wire guide Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open