



REF CV-12712

ARROW

7 Fr. | 2 Lumen | 16 cm catheter length | .032 inch dia. spring-wire guide

# Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

## Contents:

- 1: Two-Lumen Indwelling Catheter: 7 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 2: Dust Cap: Non-Vented

All components are CE 0086 unless otherwise noted.

**Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

⚠ California Prop. 65


**Warning: Cancer and Reproductive Harm**  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

Lumen	Priming Volume* (mL)		Flow Rate† (mL/hr)
	Distal (14 Ga.)	Proximal (18 Ga.)	
Distal (14 Ga.)	0.54	0.37	5300
Proximal (18 Ga.)	0.37	0.37	2100

\* Priming volumes are approximate and are done without accessories.  
† Based on 100 cm head height.



- fr** Jeu de cathétérisme veineux central à deux lumières avec cathéter Blue FlexTip
- de** Besteck für zweilumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- it** Set per cateterismo della vena centrale a due lumi con catetere Blue FlexTip
- pl** Dwukanałowy zestaw do cewnikowania żył centralnych z cewnikiem Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip
- ru** Комплект для катетеризации главных вен с двухпросветным катетером Blue FlexTip
- sl** Komplet za centralno venosko katetizacijo z dvema svetlinama s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central de dos luces con catéter Blue FlexTip
- sv** Set för tvåkanalig central venkatetrisering med Blue FlexTip kateter
- tr** İki Lümenli Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile

CE 0086  
 Teleflex Medical  
 IDA Business and Technology Park  
 Dublin Road, Athlone  
 Co. Westmeath, Ireland

Manufacturer:  
 Arrow International, Inc.  
 Subsidiary of Teleflex Incorporated  
 2400 Bernville Road  
 Reading, PA 19605 USA  
**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

SZ-12712-100C (6/18)



STERILE EO



LOT



REF CV-12712

7 Fr. | 2 Lumen | 16 cm catheter length | .032 inch dia. spring-wire guide

ARROW

Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter