



REF CV-12853



8.5 Fr. | 3 Lumen | 16 cm catheter length | .032 inch dia. spring-wire guide

# Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

## Contents:


- 1: Multi-Lumen Indwelling Catheter: 8.5 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 3: Dust Cap: Non-Vented

All components are CE 0086 unless otherwise noted.

**Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**California Prop. 65**  
**Warning: Cancer and Reproductive Harm**  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

Not made with natural rubber latex.  
Fluid path components are non-pyrogenic.




Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.41	3785
Medial (14 Ga.)	0.81	8480
Proximal (16 Ga.)	0.44	3640

\* Priming volumes are approximate and are done without accessories.  
† Based on 100 cm head height.

- cs** Vicelumenový centrální žilní katetizační set s katetrem Blue FlexTip
- da** Multilumen centralvenekateteriseringssæt med Blue FlexTip kateter
- nl** Multilumen centraalveneuze katheterisatieset met Blue FlexTip katheter
- fr** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- de** Besteck für viellumigen zentralen Venenkatheter mit Blue FlexTip Katheter
- el** Σετ καθητέρισμε veineux central à lumières multiples με καθητέρα Blue FlexTip
- it** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- no** Multilumen sentral venekateteriseringssett med Blue FlexTip kateter
- pl** Zestaw do cewnikowania żyły centralnej z cewnikiem wielokanałowym z końcówką Blue FlexTip
- pt** Conjunto de cateterização venosa central com lúmen múltiplo e cateter Blue FlexTip
- ru** Набор с многопросветным центральным венозным катетером с кончиком Blue FlexTip
- sk** Multilúmenová centrálna venózna katetizačná sada s katétrom Blue FlexTip
- sl** Večlumenski centralni venski komplet za katetizacijo s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- sv** Set för flerkanalig central venkateterisering med Blue FlexTip kateter
- tr** Multi Lümen Santral Venöz Kateterizasyon Seti, Mavi esnek uçlu Kateter ile

Manufacturer:  
Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA



**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)  
SZ-12853-112H (10/17)



REF CV-12853

8.5 Fr. | 3 Lumen | 16 cm catheter length | .032 inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

