

Dimensions:  
Width: 6-3/4" (17.2 cm)  
Length: 11-5/16" (28.7 cm)

**REF** **EU-12853-N**

# Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

# ARROW®

**Contents:**

- 1: ..... Multi-Lumen Indwelling Catheter: 8.5 Fr x 6" (16 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: ..... Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer™
- 1: ..... Introducer Needle: 18 Ga x 2-1/2" (6.35 cm) XTW
- 1: ..... Syringe: 5 ml Luer-Slip
- 1: ..... Tissue Dilator: 9 Fr (3.0 mm) x 10.2 cm
- 1: ..... SecondSite™ Adjustable Hub: Fastener
- 1: ..... SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 0086 unless otherwise noted.

Fluid path components are non-pyrogenic.




**Arrow International, Inc.**  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**EC REP** Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

**Teleflex®**  
www.teleflex.com

- be** Многопросветный набор для центральной венозной катетеризации с катетером Blue FlexTip
- hr** Višelumenski centralni venski kateter Set s Blue FlexTip kateterom
- cs** Vícélumenový centrální žilní katetrizační set s katetrem Blue FlexTip
- da** Multilumen centralvenøs kateteriseringsset med Blue FlexTip kateter
- nl** Multi-Lumen centraal veneuze katheterisatieset met Blue FlexTip katheter
- et** Mitme valendikuga tsentraalse veenikateeteri komplekt koos Blue FlexTip kateetriga
- fi** Monilumeninen keskuslaskimokatetrisarja jossa Blue FlexTip katetri
- fr** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- de** Besteck für viellumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- el** Σετ Καθετηριασμού Κεντρικής φλέβας Πολλαπλούς Αυλούς με Καθετήρα Blue FlexTip
- hu** Többlumenű, centrális vénás katéterezési szett, Blue FlexTip katéterrel
- is** Fjöl lúmena holæðaleggsett með Blue FlexTip æðalegg
- it** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- it** Daugiakanalis centrinės venos kateterizavimo mini rinkinys su „Blue FlexTip“ kateteriu
- no** Multilumen sentralt venekateteriseringssett med Blue FlexTip kateter
- pl** Zestaw do cewnikowania żyły centralnej z cewnikiem wielokanalowym z końcówką Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip
- ro** Set de cateterizare venoasă centrală multi-lumen cu cateter Blue FlexTip
- ru** Набор с многопросветным центральным венозным катетером с кончиком Blue FlexTip
- sr** Višelumenski centralni venski kateterizacioni Set sa Blue FlexTip kateterom
- sk** Viaclúmenová centrálna venózná katetrizačná súprava s katétrom Blue FlexTip
- sl** Večlumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- sv** Set för flerkanalig central venkateterisering med Blue FlexTip kateter
- tr** Multi Lumen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- uk** Багатопросвітний центральний венозний катетер Набор: катетер Blue с кінчиком Blue FlexTip Catheter



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga)	0.41	3785
Medial (14 Ga)	0.81	8480
Proximal (16 Ga)	0.44	3640

\*Priming volumes are approximate and are done without the injection cap.

†Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

SZ-12853-125B (110/16)

**LOT**



**ARROW**



**REF** EU-12853-N **8.5** Fr **3** Lumen **16** cm catheter length **.032** inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open