



REF CV-12854

ARROW®

8.5 Fr. | 4 Lumen | 16 cm catheter length | .032 inch dia. spring-wire guide

Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Contents:

- 1: Quad-Lumen Indwelling Catheter: 8.5 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 4: Dust Cap: Non-Vented

All components are CE 0086 unless otherwise noted.

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

⚠ California Prop. 65

Warning: Cancer and Reproductive Harm
www.P65Warnings.ca.gov.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

- fr** Jeu de cathétérisme veineux central à quatre lumières avec cathéter Blue FlexTip
- de** Besteck für vierlumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- it** Set per cateterismo della vena centrale a quattro lumi con catetere Blue FlexTip
- pl** Zestaw do cewnikowania żyły centralnej z cewnikiem 4-kanalowym z końcówką Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Lúmen Quádruplo e Cateter Blue FlexTip
- ru** Четырех-просветный комплект для установки центрального венозного катетера с катетером Blue FlexTip
- sl** Štiri-lumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central de cuatro luces con catéter Blue FlexTip
- sv** Set för fyrkanalig central venkatetrisering med Blue FlexTip kateter
- tr** Dört Lümen Santral Venöz Kateterizasyon Seti, Mavi esnek uçlu Kateter ile

	Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
	Distal (16 Ga.)	0.39	3200
	Medial 1 (14 Ga.)	0.47	5400
	Medial 2 (18 Ga.)	0.36	1600
	Proximal (18 Ga.)	0.38	1800

* Priming volumes are approximate and are done without accessories.
† Based on 100 cm head height.

CE 0086
Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com
SZ-12854-1021 (10/17)



LOT

Use By



REF CV-12854

8.5 Fr. | 4 Lumen | 16 cm catheter length | .032 inch dia. spring-wire guide

Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

ARROW®