

Dimensions:
Width: 6-3/4" (17,2 cm)
Length: 11-5/16" (28,7 cm)

REF EU-12854-N

Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

ARROW®

Contents:

- 1:Quad-Lumen Indwelling Catheter:
8.5 Fr x 6" (16 cm) Radiopaque
Polyurethane with Blue FlexTip®,
Extension Line Clamps
- 1: Spring-Wire Guide, Nitinol, Marked:
.032" (0.81 mm) dia. x 23-5/8"
(60 cm) (Straight Soft Tip on One
End – "J" Tip on Other) with Arrow
Advancer™ with ECG Mark
- 1: Introducer Needle: 18 Ga x 2-1/2"
(6.35 cm) XTW
- 1: Syringe: 5 ml Luer-Slip
- 1: Tissue Dilator: 9 Fr (3.0 mm)
x 10.2 cm
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter

All components are CE 0086 unless otherwise noted.


Fluid path components are non-pyrogenic.
Contains no medication.



Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

- BE** Чатырохпрасветны набор для цэнтральнай венознай катэтарызацыі з катэтрам Blue FlexTip
- HR** Quad-Lumen centralni venski kateter Komplet s Blue FlexTip kateter
- CS** Čtyřlumenový centrální žilní katetrizační set s katetrem Blue FlexTip
- DA** Firdobbelt centralvenøs kateteriseringsset med Blue FlexTip kateter
- NL** Centraal veneuze quad-lumen katheterisatieset met Blue FlexTip katheter
- ET** Neljavalendikulise tsentraalveeni kateteriseerimiskomplekt Blue FlexTip kateetriiga
- FI** Quad-Lumen keskuskaskimokatetri sarja Blue FlexTip katetrilla
- FR** Jeu de cathétérisme veineux central à quatre lumières avec cathéter Blue FlexTip
- DE** Besteck für vierlumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- EL** Σετ καθετηριασμού κεντρικής φλέβας τετραπλάσιο αυλού με καθετήρα Blue FlexTip
- HU** Quad-Lumen centrális vénás katéterezési szett, Blue FlexTip katéterrel
- IS** Fjögur lúmena holæðaleggsett með Blue FlexTip æðalegg
- IT** Set per cateterismo della vena centrale a quattro lumi con catetere Blue FlexTip
- LT** Keturių angų centrinės venos kateterizavimo rinkinys su „Blue FlexTip“ kateteriu
- NO** Sentralt venekateteriseringssett med fire lumen og Blue FlexTip kateter
- PL** Zestaw do cewnikowania żył centralnych z 4-kanalowym cewnikiem z końcówką Blue FlexTip
- PT** Conjunto de Cateterização Venosa Central com Lúmen Quádruplo e Cateter Blue FlexTip
- RO** Set de cateterizare venoasă centrală Quad-Lumen cu cateter Blue FlexTip
- RU** Четырех-просветный комплект для установки центрального венозного катетера с катетером Blue FlexTip
- SR** Quad-Lumen komplet za centralnu vensku kateterizaciju sa Blue FlexTip kateterom
- SK** Štvorlúmenová centrálna venózna katetrizačná súprava s katétrom Blue FlexTip
- SL** Štiri-lumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip
- ES** Conjunto de cateterización venosa central de cuatro luces con catéter Blue FlexTip
- SV** Set för fyrkanalyg central venkateterisering med Blue FlexTip kateter
- TR** Dörtlü-Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- UK** Чотирипросвітний комплект для встановлення центрального венозного катетера з катетером Blue FlexTip



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga)	0.39	3200
Medial 1 (14 Ga)	0.47	5400
Medial 2 (18 Ga)	0.36	1600
Proximal (18 Ga)	0.38	1800

*Priming volumes are approximate and are done without the injection cap.
†Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

LOT



(01) 0 08 01902 07880 3

57-12854-1088 (11/14)



REF EU-12854-N

8.5 Fr 4 Lumen 16 cm catheter length .032 inch dia. spring-wire guide

Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open