

Dimensions:
Width: 6-3/4" (17.2 cm)
Length: 11-5/16" (28.7 cm)

Contents:

- 1: Two-Lumen Indwelling Catheter:
4 Fr x 5-1/8" (13 cm) Radiopaque
Polyurethane with Blue FlexTip®,
Extension Line Clamps
- 1: Spring-Wire Guide, Nitinol, Marked:
.018" (0.46 mm) dia. x 17-3/4"
(45 cm) (Straight Soft Tip on One
End - "J" Tip on Other)
with Arrow Advancer™
- 1: Catheter: 22 Ga x 1-3/4"
(4.45 cm) Radiopaque over 25 Ga,
RW Introducer Needle and 5 ml
Luer-Slip Syringe
- 1: Introducer Needle: 20 Ga
x 1-1/2" (4 cm) XTW
- 1: Tissue Dilator: 5 Fr (1.8 mm)
x 7.5 cm
- 1: Tissue Dilator: 5 Fr (1.8 mm)
x 4.1 cm
- 1: SecondSite™ Adjustable Hub:
Fastener
- 1: SecondSite™ Adjustable Hub:
Catheter Clamp
- 1: Safety Scalpel: #11
- 1: ... SharpsAway II™ Locking Disposal Cup

All components are CE 0086 unless otherwise noted.

Fluid path components are non-pyrogenic.



Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA


CE REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Teleflex®
www.teleflex.com

REF EU-14402-N
Two-Lumen Central Venous Catheterization Set
with Blue FlexTip® Catheter

ARROW®

- be** Двухцветный набор для центральной венозной катетеризации с катетерам Blue FlexTip
- hr** Dvolumenski centralni venski kateter Set s Blue FlexTip kateterom
- cs** Dvoulumenový centrální žilní katetrizační set s katetrem Blue FlexTip
- da** To-lumen centralvenøs kateteriseringsset med Blue FlexTip kateter
- nl** Centraal veneuze dubbellumen katheterisatieset met Blue FlexTip katheter
- et** Kahe valendikuga tsentraalse veenikateetri komplekt koos Blue FlexTip kateetriga
- fi** Kaksilumen-keskuslaskimokatetri setti jossa Blue FlexTip katetri
- fr** Jeu de cathétérisme veineux central à deux lumières avec cathéter Blue FlexTip
- de** Besteck für zweilumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- el** Σετ Καθητηριασμού Κεντρικής φλέβας με Διπλούς Αυλούς και Καθετήρα με Blue FlexTip
- hu** Két lumenű centrális vénás katéterezési készlet Blue FlexTip katéterrel
- is** Tveggja-lúmena hölaðaleggsett með Blue FlexTip æðalegg
- it** Set per cateterismo della vena centrale a due lumi con catetere Blue FlexTip
- lt** Dviejų kanalų centrinės venos keteterizavimo rinkinys su „Blue FlexTip“ kateteriu
- no** To-Lumen Sentralt venekateteriseringssett med Blue FlexTip kateter
- pl** Zestaw do cewnikowania żył centralnych z cewnikiem 2-kanalowym z końcówką Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip
- ro** Trusa de cateterizare venoasă centrală, cu două lumene cu cateter Blue FlexTip
- ru** Двухпросветный центральный венозный катетер с наконечником Blue FlexTip
- sr** Dvo-lumenski centralno venski set za kateterizaciju sa Blue FlexTip kateterom
- sk** Dvojlúmenová katetrizačná súprava na katetrizáciu centrálnej žily s katétrom Blue FlexTip
- sl** Dvolumenski centralni venski pribor za katetrizacijo s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central de dos luces con catéter Blue FlexTip
- sv** Set för tvåkanalig central venkatetrering med Blue FlexTip kateter
- tr** Blue FlexTip Kateterli İki Lümenli Santral Venöz Kateterizasyon Seti
- uk** Двохпросвітний центральный венозный катетер з наконечником Blue FlexTip



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (22 Ga)	0.25	550
Proximal (22 Ga)	0.31	425

* Priming volumes are approximate and are done without the injection cap.
† Based on 100 cm head height.

SZ-14402-110B (9/16)

LOT



REF EU-14402-N 4 Fr 2 Lumen 13 cm catheter length .018 inch dia. spring-wire guide
Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Read to Order