



Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter



REF CV-14702-BZ

- pt** **Conjunto de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip**
Manter afastado da luz solar. Manter seco.
- es** **Conjunto de cateterización venosa central de dos lúmenes con catéter Blue FlexTip**
Mantener alejado de la luz del sol. Mantener en un lugar seco.

Contents:

- 1: Two-Lumen Indwelling Catheter: 7 Fr. x 30 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™ with ECG Mark
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator
- 1: Injection Site Cap
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 0086 unless otherwise noted.

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

Sterile. Fluid path components are non-pyrogenic.



Caution



Consult instructions for use



Do not reuse



Do not resterilize

STERILE EO

Sterilized by ethylene oxide



Keep away from sunlight



Keep dry



Do not use if package is damaged




Not made with natural rubber latex

Manufacturer
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

CE 0086 Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Teleflex®
www.teleflex.com



Catheter Cross Section

Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.57	2310
Proximal (16 Ga.)	0.62	2410

* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.

52-14702-123A (3/17)



Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

REF CV-14702-BZ

7 Fr.

2 Lumen

30 cm catheter length

.032 inch dia. spring-wire guide

