

Dimensions:
 Width: 6-3/4" (17.2 cm)
 Length: 11-5/16" (28.7 cm)

PMS 259

REF **BR-14703-EK**

Multi-Lumen Central
 Venous Catheterization Set
 with Blue FlexTip® Catheter

ARROW®

Contents:

- 1: Multi-Lumen Indwelling Catheter:
 7 Fr x 12" (30 cm) Radiopaque
 Polyurethane with Blue FlexTip®,
 Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032"
 (0.81 mm) dia. x 26-3/4" (68 cm)
 (Straight Soft Tip on One End –
 "J" Tip on Other) with Arrow
 Advancer™ with ECG Mark
- 1: Introducer Needle: 18 Ga
 x 2-1/2" (6.35 cm) XTW
 and 5 ml Luer-Slip Syringe
- 1: Tissue Dilator: 8.5 Fr (2.8 mm)
 x 10.2 cm
- 1: SecondSite™ Adjustable Hub:
 Fastener
- 1: SecondSite™ Adjustable Hub:
 Catheter Clamp
- 1: Scalpel: #11

All components are CE 0086 unless otherwise noted.

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.


Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

Sterile. Fluid path components are non-pyrogenic.

Manufacturer
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA

CE 0086
 Teleflex Medical
 IDA Business and Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland

Teleflex®
 www.teleflex.com



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga)	0.49	2300
Medial (18 Ga)	0.40	1000
Proximal (18 Ga)	0.44	1100

*Priming volumes are approximate and are done without the injection cap.
 †Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

- fr** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- de** Besteck für viellumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- it** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip
- es** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- sv** Set för flerkanalig central venkatetrisering med Blue FlexTip kateter

PMS 259

SZ14703114F (5/16)

- Caution
- Consult instructions for use
- Do not reuse
- Do not resterilize
- Sterilized by ethylene oxide
- Keep away from sunlight
- Keep dry
- Do not use if package is damaged
- Not made with natural rubber latex

LOT

Use By

PMS 259

PMS 259

ARROW®
 INTERNATIONAL



REF **BR-14703-EK**

7 Fr **3 Lumen** **30 cm** catheter length **.032** inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open