

Dimensions:  
Width: 6-3/4" (17,2 cm)  
Length: 11-5/16" (28,7 cm)

# REF EU-14703-EN

## Multi-Lumen Central Venous Catheterization Set

### with Blue FlexTip® Catheter

# ARROW®

#### Contents:

- 1: ..... Multi-Lumen Indwelling Catheter:  
7 Fr x 12" (30 cm) Radiopaque  
Polyurethane with Blue FlexTip®,  
Extension Line Clamps,
- 1: ..... Spring-Wire Guide, Nitinol,  
Marked: .032" (0.81 mm) dia.  
x 26-3/4" (68 cm) (Straight Soft  
Tip on One End - "J" Tip on Other)  
with Arrow Advancer™  
with ECG mark
- 1: ..... Introducer Needle: 18 Ga  
x 2-1/2" (6.35 cm) XTW
- 1: ..... Arrow® Raulerson Spring-Wire  
Introduction Syringe: 5 ml
- 1: ..... Pressure Transduction Probe
- 1: ..... Fastener: Catheter Clamp
- 1: ..... Clamp: Catheter
- 1: ..... Tissue Dilator 8.5 Fr (2.8 mm)  
x 4" (10.2 cm)

All components are CE 0086 unless otherwise noted.


Fluid path components are non-pyrogenic.  
Contains no medication.



 Arrow International, Inc.  
2400 Bernville Road  
Reading, PA 19605 USA

 Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

- CS** Vícelenový centrální žilní katetizační set s katetrem Blue FlexTip
- DA** Multilumen centralvenøs kateteriseringssæt med Blue FlexTip-kateter
- NL** Multi-lumen centraal veneuze katheterisatieset met Blue FlexTip katheter
- FR** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- DE** Besteck für viellumigen, zentralen Venenkatheter mit Blue-FlexTip-Katheter
- EL** Σετ Καθετηριασμού Κεντρικής φλέβας Πολλαπλούς Αυλούς με Καθετήρα Blue FlexTip
- IT** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- NO** Multilumen sentralt venekateteriseringssett med Blue FlexTip kateter
- PL** Zestaw do cewnikowania żyły centralnej z cewnikiem wielokanałowym z końcówką Blue FlexTip
- PT** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip
- RU** Набор с многопросветным центральным венозным катетером с кончиком Blue FlexTip
- SL** Večlumenski centralni venski komplet za katetizacijo s katetrom Blue FlexTip
- ES** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- SV** Set för flerkanalig central venkateterisering med Blue FlexTip-kateter
- TR** Multi Lúmen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)
Distal (16 Ga)	0.49	2300
Medial (18 Ga)	0.40	1000
Proximal (18 Ga)	0.44	1100

\* Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

5214703-137B (6/13)

LOT



(01) 0 08 01902 07934 3



REF EU-14703-EN 7 Fr 3 Lumen 30 cm catheter length .032 inch dia. spring-wire guide  
Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open