



REF MP-14703-CA

# Three-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

**7** Fr. ( **2.4** mm ) | **3** Lumen | **30** cm catheter length | **.018** inch dia. spring-wire guide


**Contents:**

- |  |   |
|--|---|
| 1: Three-Lumen Indwelling Catheter: 7 Fr. x 12" (30 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps and Stylet                                 | 1: Syringe: 2 mL Luer-Slip                    |
| 1: Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 27-1/2" (70 cm) (Angled Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer and hub | 1: Tissue Dilator: 8 Fr. x 9.7 cm             |
| 1: Catheter: 20 Ga. (0.84 mm) x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. (0.72 mm) TW Introducer Needle   | 3: Dust Cap: Non-Vented                       |
| 1: Introducer Needle: Echogenic 22 Ga. (0.72 mm) x 2-3/4" (7 cm) XTW   | 1: SecondSite™ Adjustable Hub: Fastener       |
| 1: Introducer Needle: Echogenic 22 Ga. (0.72 mm) x 1-3/8" (3.5 cm) XTW   | 1: SecondSite™ Adjustable Hub: Catheter Clamp |
|  | 2: Guidewire Introducer                       |
|  | 1: Y-hub                                      |
|  | 3: Towel: 13" x 18" (33 cm x 46 cm)           |

**Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.49	2300
Medial (18 Ga.)	0.40	1000
Proximal (18 Ga.)	0.44	1100

\* Priming volumes are approximate and are done without accessories.  
 † Based on 100 cm head height.



STERILE EO



LOT



Arrow International, Inc.  
 Subsidiary of Teleflex Incorporated  
 2400 Bernville Road  
 Reading, PA 19605 USA

**Teleflex®**  
 www.teleflex.com

SZ-14703-148B (09/19)



REF MP-14703-CA

**7** Fr. ( **2.4** mm ) | **3** Lumen | **30** cm catheter length | **.018** inch dia. spring-wire guide

Three-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter