


Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga)	0.56	2935
Medial (14 Ga)	1.06	5920
Proximal (16 Ga)	0.55	2515

\*Priming volume are approximate and done without the stopcock.

†Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



**Contents:**

- 1: ..... Multi-Lumen Indwelling Catheter: 8.5 Fr x 12" (30 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: ..... Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™ with ECG Mark
- 1: ..... Introducer Needle: 18 Ga x 2-1/2" (6.35 cm) XTW and 5 ml Luer-Slip Syringe
- 1: ..... Tissue Dilator: 9 Fr (3.0 mm) x 10.2 cm
- 1: ..... Fastener: Catheter Clamp
- 1: ..... Clamp: Catheter
- 1: ..... Scalpel: #11
- 3: ..... Stopcocks: 4-way High-Flow

All components are CE 0086 unless otherwise noted.

Rx only.

**Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

**Sterile.** Single use – do not resterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Not made with natural rubber latex. Keep away from sunlight. Keep dry. Fluid path components are non-pyrogenic. Contains no medication.



Arrow International, Inc.  
2400 Bernville Road  
Reading, PA 19605 USA

**EC REP** Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

**REF BR-14853-EK**

# Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

# ARROW®

- CS** Vícelumenový centrální žilní katetizační set s katetrem Blue FlexTip
- NO** Multilumen sentralt venekateteriseringsset med Blue FlexTip kateter
- DA** Multilumen centralvenekateteriseringsæt med Blue FlexTip kateter
- PL** Zestaw wielokanałowy do cewnikowania żył centralnych z cewnikiem typu Blue FlexTip
- NL** Multilumen centraal-veneuze katheterisatieset met Blue FlexTip katheter
- PT** Conjunto de cateterização venosa central com lúmen múltiplo e cateter Blue FlexTip
- FR** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- SK** Multilúmenová centrálna venózná katetizačná sada s katétrom Blue FlexTip
- DE** Besteck für viellumigen zentralen Venenkatheter mit Blue FlexTip Katheter
- ES** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- EL** Σετ κεντρικού φλεβικού καθετηριασμού πολλαπλών αυλών με καθετήρα Blue FlexTip
- SV** Set för flerkanalig central venkatetrering med Blue FlexTip kateter
- IT** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- TW** 帶 Blue FlexTip 導管的多腔中央靜脈插管套裝
- JA** マルチルーメン中心静脈カテーテルセット - Blue FlexTipカテーテル付
- TR** Multi-Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile

PMS 259

SZ-14853-104F (9/1/16)

PMS 259