

Dimensions:
Width: 6-3/4" (17,2 cm)
Length: 11-5/16" (28,7 cm)

REF EU-14854-EN

Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

ARROW®

Contents:

- 1: Quad-Lumen Indwelling Catheter: 8.5 Fr x 12" (30 cm)
Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer™ with ECG mark
- 1: Introducer Needle: 18 Ga x 2-1/2" (6.35 cm) XTW
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 ml
- 1: Pressure Transduction Probe
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter
- 1: Tissue Dilator 10 Fr (3.5 mm) x 4" (10.2 cm)

All components are CE 0086 unless otherwise noted.

Fluid path components are non-pyrogenic. Contains no medication.



Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

CS Čtyřlumenový centrální žilní katetrizační set s katetrem Blue FlexTip

DA Firdobbelt centralvenøs kateteriseringssæt med Blue FlexTip kateter

NL Centraal veneuze vier-lumen katheterisatieset met Blue FlexTip katheter

FR Jeu de cathétérisme veineux central à quatre lumières avec cathéter Blue FlexTip

DE Besteck für vierlumigen, zentralen Venenkatheter mit Blue-FlexTip-Katheter

EL Σετ καθετηριασμού κεντρικής φλέβας τετραπλάσιο αυλού με καθετήρα Blue FlexTip

IT Set per cateterismo della vena centrale a quattro lumi con catetere Blue FlexTip

NO Sentralt venekateteriseringssett med fire lumen og Blue FlexTip kateter

PL Zestaw do cewnikowania żył centralnych z czterokanałowym cewnikiem z końcówką Blue FlexTip

PT Conjunto de Cateterização Venosa Central com Lúmen Quádruplo e Cateter Blue FlexTip


RU Четырех-просветный комплект для установки центрального венозного катетера с катетером Blue FlexTip

SL Štiri-lumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip

ES Conjunto de cateterización venosa central de cuatro lúmenes con catéter Blue FlexTip

SV Set för fyrkanalyg central venkatetrering med Blue FlexTip-kateter

TR Dörtlü-Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)
Distal (16 Ga)	0.52	2100
Medial 1 (14 Ga)	0.66	3900
Medial 2 (18 Ga)	0.44	1000
Proximal (18 Ga)	0.47	1000

*Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.

†Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

LOT



(01) 0 08 01902 07936 7

52-14854-1048 (6/13)



REF EU-14854-EN 8.5 Fr 4 Lumen 30 cm catheter length .032 inch dia. spring-wire guide
Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open