

Dimensions: 8-5/8" (21,9 cm) x 15-1/2" (39,4 cm)

REF EU-14854-CVT

**Quad-Lumen
Central Venous Catheterization Set
with Blue FlexTip® Catheter**

ARROW®

Contents:

- 1: Quad-Lumen Indwelling Catheter: 8.5 Fr x 12" (30 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™ with ECG Mark
- 1: Introducer Needle: Echogenic 18 Ga x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 ml Luer-Slip
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter
- 1: Tissue Dilator 9 Fr (3.0 mm) x 4" (10.2 cm)
- 1: ECG Extension Cable
- 1: Scalpel: #11
- 3: Injection Site Cap

All components are CE 0086 unless otherwise noted.


Fluid path components are non-pyrogenic. Contains no medication.



Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

- CS** Čtyřlumenový centrální žilní katetrizační set s katetrem Blue FlexTip
- DA** Firdobbelt centralvenos kateteriseringsset med Blue FlexTip kateter
- NL** Centraal veneuze vier-lumen katheterisatieset met Blue FlexTip katheter
- ET** Neljavalendikulise tsentraalveeni kateteriseerimiskomplekt Blue FlexTip katetriga
- FI** Quad-Lumen keskuslaskimokatetrisarja Blue FlexTip katetrilla
- FR** Jeu de cathétérisme veineux central à quatre lumières avec cathéter Blue FlexTip
- DE** Besteck für vierlumigen, zentralen Venenkatheter mit Blue-FlexTip-Katheter
- EL** Σετ καθετηριασμού κεντρικής φλέβας τετραπλάσιο αυλού με καθετήρα Blue FlexTip
- HU** Quad-Lumen centrális vénás katéterezési szett, Blue FlexTip katéterrel
- IT** Set per cateterismo della vena centrale a quattro lumi con catetere Blue FlexTip
- NO** Sentralt venekateteriseringssett med fire lumen og Blue FlexTip kateter
- PL** Zestaw do cewnikowania żył centralnych, z czterokanałowym cewnikiem z końcówką Blue FlexTip
- PT** Conjunto de Cateterização Venosa Central com Lúmen Quádruplo e Cateter Blue FlexTip
- SK** Štvorlúmenová centrálna venózna katetrizačná súprava s katétrom Blue FlexTip
- ES** Conjunto de cateterización venosa central de cuatro luces con catéter Blue FlexTip
- SV** Set för fyrkanalyg central venkatetrering med Blue FlexTip-kateter



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)
Distal (16 Ga)	0.52	2100
Medial 1 (14 Ga)	0.66	3900
Medial 2 (18 Ga)	0.44	1000
Proximal (18 Ga)	0.47	1000

*Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.

†Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

LOT



(01) 0 08 01902 08019 6

SE14854107A(5/13)



REF EU-14854-CVT **8.5 Fr** 4 Lumen **30 cm** catheter length **.032** inch dia. spring-wire guide
Quad-Lumen Central Venous Catheterization Set with Blue® FlexTip Catheter

Peel to Open