

Dimensions:  
Width: 6-3/4" (17.2 cm)  
Length: 11-5/16" (28.7 cm)

**REF EU-15403-N**  
**Multi-Lumen Central Venous Catheterization Set**  
with Blue FlexTip® Catheter  
**ARROW®**

**Contents:**

- 1: ..... Multi-Lumen Indwelling Catheter: 4 Fr x 3" (8 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: ..... Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer™
- 1: ..... Catheter: 22 Ga x 1-3/4" (4.45 cm) Radiopaque over 25 Ga, RW Introducer Needle and 5 ml Luer-Slip Syringe
- 1: ..... Introducer Needle: 20 Ga x 1-1/2" (4 cm) XTW
- 1: ..... Tissue Dilator 5 Fr (1.8 mm) x 3" (7.5 cm)
- 1: ..... Tissue Dilator 5 Fr (1.8 mm) x 1-5/8" (4.1 cm)
- 1: ..... Fastener: Catheter Clamp
- 1: ..... Clamp: Catheter
- 1: ..... Safety Scalpel: #11
- 1: ... SharpsAway II™ Locking Disposal Cup

All components are CE 0086 unless otherwise noted.


Contains no medication.  
Fluid path components are non-pyrogenic.



Arrow International, Inc.  
2400 Bernville Road  
Reading, PA 19605 USA

**EC/REP** Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

- BE** Многопросветный набор для центральной венозной катетеризации з катэтрам Blue FlexTip
- HR** Višelumenski centralni venski kateter Set s Blue FlexTip kateterom
- CS** Vícélumenový centrální žilní katetrizační set s katetrem Blue FlexTip
- DA** Multilumen centralvenøs kateteriseringssæt med Blue FlexTip kateter
- NL** Multi-lumen centraal veneuze katheterisatieset met Blue FlexTip katheter
- ET** Mitme valendikuga tsentraalse veenikateeteri komplekt koos Blue FlexTip kateetriga
- FI** Monilumeninen keskuslaskimokatetri sarja jossa sininen FlexTip katetri
- FR** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- DE** Besteck für viellumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- EL** Σετ Καθετηριασμού Κεντρικής φλέβας Πολλαπλούς Αυλούς με Καθετήρα Blue FlexTip
- HU** Többlumenű, centrális vénás katéterezési szett, Blue FlexTip katéterrel
- IS** Fjöllúmena hölaðaleggsett með Blue FlexTip æðalegg
- IT** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- LT** Daugiakanaliscentrinis venos kateterizavimo mini rinkinys su „Blue FlexTip“ kateteriu
- NO** Multilumen sentralt venekateriseringssett med Blue FlexTip kateter
- PL** Zestaw do cewnikowania żyły centralnej z cewnikiem wielokanalowym z końcówką Blue FlexTip
- PT** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip
- RO** Set de cateterizare venoasă centrală multi-lumen cu cateter Blue FlexTip
- RU** Набор с многопросветным центральным венозным катетером с кончиком Blue FlexTip
- SR** Višelumenski centralni venski kateterizacioni Set sa Blue FlexTip kateterom
- SK** Viaclúmenová centrálna venózna katetrizačná súprava s katétrom Blue FlexTip
- SL** Večlumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip
- ES** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- SV** Set för flerkanalig central venkateterisering med Blue FlexTip kateter
- TR** Multi Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- UK** Багатопросвітний центральний венозний катетер Набор: катетер Blue с кінчиком Blue FlexTip Catheter



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)
Distal (20 Ga)	0.27	1020
Medial (23 Ga)	0.29	255
Proximal (23 Ga)	0.32	300

\* Priming volumes are approximate and are done without the injection cap.  
† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

SZ-15403-102A (2/14)



(01) 0 08 01902 08773 7



**REF EU-15403-N**

4 Fr 3 Lumen

8 cm catheter length

.018 inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open