

Dimensions:  
Width: 10-5/8" (27 cm)  
Length: 12-1/2" (31,8 cm)

# REF AK-15703-J

## Multi-Lumen Central Venous Catheterization Kit


with Blue FlexTip® Catheter

# ARROW®

Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)
Distal (16 Ga)	0.44	3100
Medial (18 Ga)	0.39	1500
Proximal (18 Ga)	0.39	1600

\* Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



Fluid path components are non-pyrogenic. Contains no medication.



Arrow International, Inc.  
2400 Bernville Road, Reading, PA 19605 USA

EC REP Teleflex Medical  
IDA Business and Technology Park  
Athlone, Ireland

### Special Order Product

#### Contents:

- 1: Multi-Lumen Indwelling Catheter: 7 Fr x 8" (20 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps, Injection Site Caps
- 1: Spring-Wire Guide, Marked: .032" (.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™
- 1: Catheter: 18 Ga x 2-1/2" (6.35 cm) Radiopaque over 20 Ga RW Introducer Needle
- 1: Fastener: Catheter Clamp
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 22 Ga x 2-3/4" (6.98 cm) TW
- 1: Introducer Needle: 18 Ga x 2-1/2" (6.35 cm) XTW and 5 ml Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: 22 Ga x 1-1/2" (3.81 cm) and 5 ml Luer-Slip Syringe
- 1: Injection Needle: 25 Ga x 1" (2.54 cm) and 3 ml Luer-Lock Syringe
- 3: Dust Caps: Non-Vented
- 1: Clamp: Catheter
- 1: Tissue Dilator
- 1: CSR Wrap
- 1: Drape: 21" x 36" (53 x 91 cm) with 3" (7.6 cm) fenestration, with adhesive
- 5: Gauze Pads: 4" x 4" (10 x 10 cm)
- 2: Gauze Pads: 2" x 2" (5 x 5 cm)
- 1: Scalpel: #11
- 2: Prep Sponge Swabs
- 3: Stopcocks: 3-way
- 1: Suture: 3-0 Silk with Straight Needle
- 1: Tray: Prep
- 3: UserGuard® Needle-Free Injection Hubs
- 1: SharpsAway® Disposal Cup

SZ157031180 (10/11)

REF **AK-15703-J** 7 Fr 3 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Kit with Blue FlexTip®



(01) 0 08 01902 00036 1

**ARROW**  
INTERNATIONAL  
Special Order Product

REF **AK-15703-J** 7 Fr 3 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Kit with Blue FlexTip®

Peel to Open