

Dimensions:
Width: 6-3/4" (17.2 cm)
Length: 11-5/16" (28.7 cm)

REF **EU-15703-CVT**
Multi-Lumen Central Venous Catheterization Set
with Blue FlexTip® Catheter
ARROW®

Contents:

- 1: Multi-Lumen Indwelling Catheter:
7 Fr x 8" (20 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™ with ECG Mark
- 1: Introducer Needle: Echogenic 18 Ga x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 ml Luer-Slip
- 1: Tissue Dilator: 8.5 Fr (2.8 mm) x 10.2 cm
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter
- 1: Scalpel: #11
- 1: ECG Extension Cable
- 2: Injection Site Cap

All components are CE 0086 unless otherwise noted.

Contains no medication.


Fluid path components are non-pyrogenic.



Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

- BE** Многопросветный набор для центральной венозной катетеризации з катетэрам Blue FlexTip
- HR** Višelumenski centralni venski kateter Set s Blue FlexTip kateterom
- CS** Vícelumenový centrální žilní katetrizační set s katetrem Blue FlexTip
- DA** Multilumen centralvenøs kateteriseringsset med Blue FlexTip kateter
- NL** Multi-lumen centraal veneuze katheterisatieset met Blue FlexTip katheter
- ET** Mitme valendikuga tsentraalse veenikateeteri komplekt koos Blue FlexTip kateetriga
- F** Moniluumeninen keskuslaskimokatetri sarja jossa sininen FlexTip katetri
- FR** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- DE** Besteck für viellumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- EL** Σετ Καθετηριασμού Κεντρικής φλέβας Πολλαπλούς Αυλούς με Καθετήρα Blue FlexTip
- HU** Többlumenű, centrális vénás katéterezési szett, Blue FlexTip katéterrel
- IS** Fjöllumena holæðaleggsett með Blue FlexTip æðalegg
- IT** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- LT** Daugiakanaliscentrinës venos kateterizavimo mini rinkinys su „Blue FlexTip“ kateteriu
- NO** Multilumen sentralt venekateriseringssett med Blue FlexTip kateter
- PL** Zestaw do cewnikowania żyły centralnej z cewnikiem wielokanałowym z końcówką Blue FlexTip
- PT** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip
- RO** Set de cateterizare venoasă centrală multi-lumen cu cateter Blue FlexTip
- RU** Набор с многопросветным центральным венозным катетером с кончиком Blue FlexTip
- SR** Višelumenski centralni venski kateterizacioni Set sa Blue FlexTip kateterom
- SK** Vialúmenová centrálna venózna katetrizačná súprava s katétrom Blue FlexTip
- SL** Večlumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip
- ES** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- SV** Set för flerkanalig central venkatetrisering med Blue FlexTip kateter
- TR** Multi Lúmen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- UK** Багатопроевітний центральний венозний катетер Набір: катетер Blue с кінчиком Blue FlexTip Catheter



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga)	0.44	3100
Medial (18 Ga)	0.39	1500
Proximal (18 Ga)	0.39	1600

* Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

SZ-15703-1248 (11/14)

LOT



(01) 0 08 01902 06576 6

ARROW
INTERNATIONAL



REF **EU-15703-CVT**

7 Fr 3 Lumen

20 cm catheter length

.032 inch dia. spring-wire guide

Peel to Open

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter