

Dimensions:
Width: 6-3/4" (17.2 cm)
Length: 11-5/16" (28.7 cm)

REF **EU-15703-EN**
Multi-Lumen Central Venous
Catheterization Set
with Blue FlexTip® Catheter
ARROW

Contents:

- 1: Multi-Lumen Indwelling Catheter:
7 Fr x 8" (20 cm) Radiopaque
Polyurethane with Blue FlexTip®,
Extension Line Clamps
- 1: Spring-Wire Guide, Nitinol,
Marked: .032" (0.81 mm)
dia. x 23-5/8" (60 cm) (Straight
Soft Tip on One End – "J" Tip
on Other) with Arrow Advancer™
with ECG Mark
- 1: Introducer Needle: 18 Ga
x 2-1/2" (6.35 cm) XTW
- 1: Arrow® Raulerson Spring-Wire
Introduction Syringe: 5 ml
- 1: Tissue Dilator: 8.5 Fr (2.8 mm)
x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter
Clamp
- 1: Pressure Transduction Probe

All components are CE 0086 unless otherwise noted.

Fluid path components are non-pyrogenic.




Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Teleflex
www.teleflex.com

- be** Многопросветный набор для центральной венозной катетеризации с катетером Blue FlexTip
- hr** Višelumenski centralni venski kateter Set s Blue FlexTip kateterom
- cs** Vícélumenový centrální žilní katetrizační set s katetrem Blue FlexTip
- da** Multilumen centralvenøs kateteriseringsset med Blue FlexTip kateter
- nl** Multi-Lumen centraal veneuze katheterisatieset met Blue FlexTip katheter
- et** Mitme valendikuga tsentraalse veenikateeteri komplekt koos Blue FlexTip kateetriga
- fi** Monilumeninen keskuslaskimokatetrisarja jossa Blue FlexTip katetri
- fr** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- de** Besteck für viellumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- el** Σετ Καθητηριασμού Κεντρικής φλέβας Πολλαπλούς Αυλούς με Καθέτηρα Blue FlexTip
- hu** Többlumenű, centrális vénás katéterezési szett, Blue FlexTip katéterrel
- is** Fjöl lúmena holæðaleggsett með Blue FlexTip æðalegg
- it** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- it** Daugiakanalis centrinės venos kateterizavimo mini rinkinys su „Blue FlexTip“ kateteriu
- no** Multilumen sentralt venekateteriseringssett med Blue FlexTip kateter
- pl** Zestaw do cewnikowania żyły centralnej z cewnikami wielokanałowymi z końcówką Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip
- ro** Set de cateterizare venoasă centrală multi-lumen cu cateter Blue FlexTip
- ru** Набор с многопросветным центральным венозным катетером с кончиком Blue FlexTip
- sr** Višelumenski centralni venski kateterizacioni Set sa Blue FlexTip kateterom
- sk** Vialúmenová centrálna venózna katetrizačná súprava s katétrom Blue FlexTip
- sl** Večlumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- sv** Set för flerkanalig central venkateterisering med Blue FlexTip kateter
- tr** Multi Lumen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- uk** Багатопроевітний центральний венозний катетер Набор: катетер Blue с кінчиком Blue FlexTip Catheter



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga)	0.44	3100
Medial (18 Ga)	0.39	1500
Proximal (18 Ga)	0.39	1600

*Priming volumes are approximate and are done without the injection cap.

†Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

SZ-15703-1368 (110/116)

LOT



ARROW



REF EU-15703-EN **7** Fr **3** Lumen **20** cm catheter length **.032** inch dia. spring-wire guide
Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open