

Dimensions:  
Width: 6-3/4" (17.2 cm)  
Length: 11-5/16" (28.7 cm)

**REF EU-15703-N**  
Multi-Lumen Central Venous  
Catheterization Set  
with Blue FlexTip® Catheter

# ARROW®

**Contents:**

- 1: ..... Multi-Lumen Indwelling Catheter:  
7 Fr x 8" (20 cm) Radiopaque  
Polyurethane with Blue FlexTip®,  
Extension Line Clamps
- 1: ..... Spring-Wire Guide, Nitinol,  
Marked: .032" (0.81 mm)  
dia. x 23-5/8" (60 cm) (Straight  
Soft Tip on One End – "J" Tip  
on Other) with Arrow Advancer™  
with ECG Mark
- 1: ..... Introducer Needle: 18 Ga  
x 2-1/2" (6.35 cm) XTW
- 1: ..... Syringe: 5 ml Luer-Slip
- 1: ..... Tissue Dilator: 8.5 Fr (2.8 mm)  
x 10.2 cm
- 1: ..... SecondSite™ Adjustable Hub: Fastener
- 1: ..... SecondSite™ Adjustable Hub:  
Catheter Clamp

All components are CE 0086 unless otherwise noted.

Fluid path components are non-pyrogenic.



Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA


**EC REP** Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

**Teleflex®**  
www.teleflex.com

**LOT**



- be** Многопросветный набор для центральной венозной катетеризации с катетром Blue FlexTip
- hr** Višelumenski centralni venski kateter Set s Blue FlexTip katetrom
- cs** Vícélumenový centrální žilní katetrizační set s katetrem Blue FlexTip
- da** Multilumen centralvenøs kateteriseringsset med Blue FlexTip kateter
- nl** Multi-Lumen centraal veneuze katheterisatieset met Blue FlexTip katheter
- et** Mitme valendikuga tsentraalse veenikateeteri komplekt koos Blue FlexTip kateetriga
- fi** Monilumeninen keskuslaskimokatetrisarja jossa Blue FlexTip katetri
- fr** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- de** Besteck für viellumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- el** Σετ Καθετηριασμού Κεντρικής φλέβας Πολλαπλούς Αυλούς με Καθετήρα Blue FlexTip
- hu** Többlumenű, centrális vénás katéterezési szett, Blue FlexTip katéterrel
- is** Fjöl lúmena holæðaleggsett með Blue FlexTip æðalegg
- it** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- it** Daugiakanalis centrinės venos kateterizavimo mini rinkinys su „Blue FlexTip“ kateteriu
- no** Multilumen sentralt venekateriseringssett med Blue FlexTip kateter
- pl** Zestaw do cewnikowania żyły centralnej z cewnikiem wielokanałowym z końcówką Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip
- ro** Set de cateterizare venoasă centrală multi-lumen cu cateter Blue FlexTip
- ru** Набор с многопросветным центральным венозным катетером с кончиком Blue FlexTip
- sr** Višelumenski centralni venski kateterizacioni Set sa Blue FlexTip katetrom
- sk** Viaclúmenová centrálna venózna katetrizačná súprava s katétrom Blue FlexTip
- sl** Večlumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- sv** Set för flerkanalig central venkateterisering med Blue FlexTip kateter
- tr** Multi Lúmen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- uk** Багатопросвітний центральний венозний катетер Набор: катетер Blue с кінчиком Blue FlexTip Catheter



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga)	0.44	3100
Medial (18 Ga)	0.39	1500
Proximal (18 Ga)	0.39	1600

\*Priming volumes are approximate and are done without the injection cap.

†Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

52-15703-1378 (10/16)

**ARROW**



**REF** EU-15703-N

7 Fr 3 Lumen

20 cm catheter length

.032 inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open