

Dimensions:  
Width: 8-3/4" (22.2 cm)  
Length: 11-1/8" (28.3 cm)

**REF EU-15703-CVTS**  
**Multi-Lumen Central Venous  
Catheterization Set  
with Blue FlexTip® Catheter**

# ARROW®

**Contents:**

- 1: ..... Multi-Lumen Indwelling Catheter:  
7 Fr x 8" (20 cm) Radiopaque  
Polyurethane with Blue FlexTip®,  
Extension Line Clamps
- 1: ..... Spring-Wire Guide, Marked: .032"  
(0.81 mm) dia. x 23-5/8" (60 cm)  
(Straight Soft Tip on One End – "J" Tip  
on Other) with Arrow Advancer™  
with ECG Mark
- 1: ..... Introducer Needle: Echogenic 18 Ga  
x 2-1/2" (6.35 cm) XTW
- 1: ..... Syringe: 5 ml Luer-Slip
- 1: ..... Tissue Dilator: 8.5 Fr (2.8 mm)  
x 10.2 cm
- 1: ..... Fastener: Catheter Clamp
- 1: ..... Clamp: Catheter
- 1: ..... Safety Scalpel: #11
- 1: ..... SharpsAway II™ Locking Disposal Cup
- 1: ..... ECG Extension Cable
- 2: ..... Injection Site Cap

All components are CE 0086 unless otherwise noted.


Fluid path components are non-pyrogenic.  
Contains no medication.



Arrow International, Inc.  
2400 Bernville Road  
Reading, PA 19605 USA

EC REP Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

- BE** Многопросветный набор для центральной венозной катетеризации с катетерам Blue FlexTip
- HR** Višelumenski centralni venski kateter Set s Blue FlexTip kateterom
- CS** Vícelumenový centrální žilní katetrizační set s katetrem Blue FlexTip
- DA** Multilumen centralvenos kateteriserings-sæt med Blue FlexTip-kateter
- NL** Meerlumen centraal veneuze katheterisatieset met Blue FlexTip katheter
- ET** Mitme valendikuga tsentraalse veenikateeteri komplekt koos Blue FlexTip kateetriga
- F** Monilumeninen keskuslaskimokatetri sarja jossa Blue FlexTip-katetri
- FR** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- DE** Besteck für viellumigen, zentralen Venenkatheter mit Blue FlexTip Kathete
- EL** Σετ Καθετηριασμού Κεντρικής φλέβας Πολλαπλούς Αυλούς με Καθετήρα Blue FlexTip
- HU** Többlumenű, centrális vénás katéterezési szett, Blue FlexTip katéterrel
- IS** Fjöllúmena holæðaleggsett með Blue FlexTip æðalegg
- IT** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- LT** Daugiakanalis centrinės venos kateterizavimo mini rinkinys su Blue FlexTip kateteriu
- NO** Multilumen sentralt venekateriseringssett med Blue FlexTip kateter
- PL** Zestaw do cewnikowania żyły centralnej z cewnikiem wielokanałowym z końcówką Blue FlexTip i powłoką
- PT** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip
- RO** Set de cateterizare venoasă centrală multi-lumen cu cateter Blue FlexTip
- RU** Набор с многопросветным центральным венозным катетером с кончиком Blue FlexTip
- SR** Višelumenski centralni venski kateterizacioni Set sa Blue FlexTip kateterom
- SK** Víaclúmenová centrálna venózná katetrizačná súprava s katétrom Blue FlexTip
- SL** Večlumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip
- ES** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- SV** Set för flerkanalig central venkateterisering med Blue FlexTip-kateter
- TR** Multi Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- UK** Багатопрорітний центральний венозний катетер Набор: катетер Blue с кінчиком Blue FlexTip Catheter



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)
Distal (16 Ga)	0.44	3100
Medial (18 Ga)	0.39	1500
Proximal (18 Ga)	0.39	1600

\* Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

LOT



(01) 0 08 01902 07354 9

SL-15703-108A (12/14)



**REF EU-15703-CVTS** 7 Fr 3 Lumen 20 cm catheter length .032 inch dia. spring-wire guide  
Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open