



**REF** MP-15703-CA

# Three-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

**7** Fr. ( **2.4** mm ) | **3** Lumen | **20** cm catheter length | **.018** inch dia. spring-wire guide


**Contents:**

- 1: Three-Lumen Indwelling Catheter: 7 Fr. x 8" (20 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps and Stylet
- 1: Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 23-5/8" (60 cm) (Angled Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer and hub
- 1: Catheter: 20 Ga. (0.84 mm) x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. (0.72 mm) TW Introducer Needle
- 1: Introducer Needle: Echogenic 22 Ga. (0.72 mm) x 2-3/4" (7 cm) XTW
- 1: Introducer Needle: Echogenic 22 Ga. (0.72 mm) x 1-3/8" (3.5 cm) XTW
- 1: Syringe: 2 mL Luer-Slip
- 1: Tissue Dilator: 8 Fr. x 9.7 cm
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 2: Guidewire Introducer
- 1: Y-hub
- 3: Towel: 13" x 18" (33 cm x 46 cm)

**Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.44	3100
Medial (18 Ga.)	0.39	1500
Proximal (18 Ga.)	0.39	1600

\* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.



**LOT**



Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex®**  
www.teleflex.com

SZ-15703-153B (09/19)



**REF** MP-15703-CA

**7** Fr. ( **2.4** mm ) | **3** Lumen | **20** cm catheter length | **.018** inch dia. spring-wire guide

Three-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter