

Dimensions:
Width: 6-3/4" (17,2 cm)
Length: 11-5/16" (28,7 cm)

REF **EU-15802-CVT**
Two-Lumen Central Venous Catheterization Set
with Blue FlexTip® Catheter
ARROW®

Contents:

- 1: Two-Lumen Indwelling Catheter:
8 Fr x 8" (20 cm) Radiopaque
Polyurethane with Blue FlexTip®,
Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032"
(0.81 mm) dia. x 23-5/8" (60 cm)
(Straight Soft Tip on One End – "J" Tip
on Other) with Arrow Advancer™
with ECG Mark
- 1: Introducer Needle: Echogenic
18 Ga x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 ml Luer-Slip
- 1: Tissue Dilator: 9 Fr (3.0 mm)
x 10.2 cm
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter
- 1: Scalpel: #11
- 1: ECG Extension Cable
- 2: Injection Site Cap

All components are CE 0086 unless otherwise noted.


Fluid path components are non-pyrogenic.
Contains no medication.



Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

- BE** Двухпросветный набор для центральной венозной катетеризации с катетером Blue FlexTip
- HR** Dvolumenski centralni venski kateter Set s Blue FlexTip kateterom
- CS** Dvoulumenový centrální žilní katetrizační set s katetrem Blue FlexTip
- DA** To-lumen centralvenøs kateteriseringsset med Blue FlexTip kateter
- NL** Centraal veneuze dubbellumen katheterisatieset met Blue FlexTip katheter
- ET** Kahe valendikuga tsentraalse veenikateetri komplekt koos Blue FlexTip kateetriga
- FI** Kaksiluumen-keskuslaskimokatetri setti jossa Blue FlexTip katetri
- FR** Jeu de cathétérisme veineux central à deux lumières avec cathéter Blue FlexTip
- DE** Besteck für zweilumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- EL** Σετ Καθητηριασμού Κεντρικής φλέβας με Διπλούς Αυλούς και Καθετήρα με Blue FlexTip
- HU** Két lumenű centrális vénás katéterezési készlet Blue FlexTip katéterrel
- IS** Tveggja-lúmena holæðaleggsett með Blue FlexTip æðalegg
- IT** Set per cateterismo della vena centrale a due lumi con catetere Blue FlexTip
- LT** Dviejų kanalų centrinės venos keteterizavimo rinkinys su „Blue FlexTip“ kateteriu
- NO** To-Lumen Sentralt venekateteriseringssett med Blue FlexTip kateter
- PL** Zestaw do cewnikowania żył centralnych z cewnikiem 2-kanalowym z końcówką Blue FlexTip
- PT** Conjunto de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip
- RO** Trusa de cateterizare venoasă centrală, cu două lumene cu cateter Blue FlexTip
- RU** Двухпросветный центральный венозный катетер с наконечником Blue FlexTip
- SR** Dvo-lumenski centralno venski set za kateterizaciju sa Blue FlexTip kateterom
- SK** Dvojlúmenová katetrizačná súprava na katetrizáciu centrálnej žily s katétrom Blue FlexTip
- SL** Dvolumenski centralni venski pribor za katetrizacijo s katetrom Blue FlexTip
- ES** Conjunto de cateterización venosa central de dos luces con catéter Blue FlexTip
- SV** Set för tvåkanalig central venkatetrisering med Blue FlexTip kateter
- TR** Blue FlexTip Kateterli İki Lümenli Santral Venöz Kateterizasyon Seti
- UK** Двохпросвітний центральный венозный катетер з наконечником Blue FlexTip



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (14 Ga)	0.75	5600
Proximal (14 Ga)	0.80	5100


* Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.
† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

5Z-15802-1128 (10/14)

LOT



(01) 0 08 01902 06581 0



REF

EU-15802-CVT

8 Fr

2 Lumen

20 cm
catheter length

.032
inch dia.

spring-wire guide

Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open