

Dimensions:  
Width: 6-3/4" (17,2 cm)  
Length: 11-5/16" (28,7 cm)

**REF EU-15853-N**  
**Multi-Lumen Central Venous Catheterization Set**  
with Blue FlexTip® Catheter

# ARROW®

**Contents:**

- 1: ..... Multi-Lumen Indwelling Catheter: 8.5 Fr x 8" (20 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: ..... Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™ with ECG Mark
- 1: ..... Introducer Needle: 18 Ga x 2-1/2" (6.35 cm) XTW and 5 ml Luer-Slip Syringe
- 1: ..... Tissue Dilator: 9 Fr (3.0 mm) x 10.2 cm
- 1: ..... Fastener: Catheter Clamp
- 1: ..... Clamp: Catheter

All components are CE 0086 unless otherwise noted.


Fluid path components are non-pyrogenic. Contains no medication.



Arrow International, Inc.  
2400 Bemville Road  
Reading, PA 19605 USA

**EC REP** Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

- BE** Многосветный набор для центральной венозной катетеризации с катетером Blue FlexTip
- HR** Višelumenski centralni venski kateter Set s Blue FlexTip kateterom
- CS** Vícélumenový centrální žilní katetrizační set s katetrem Blue FlexTip
- DA** Multilumen centralvenøs kateteriseringsset med Blue FlexTip kateter
- NL** Multi-lumen centraal veneuze katheterisatieset met Blue FlexTip katheter
- ET** Mitme valendikuga tsentraalse veenikateeteri komplekt koos Blue FlexTip kateetriga
- FI** Moniluumeninen keskuskaskimokatetri sarja jossa sininen FlexTip katetri
- FR** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- DE** Besteck für viellumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- EL** Σετ Καθητηριασμού Κεντρικής φλέβας Πολλαπλούς Αυλούς με Καθήτηρα Blue FlexTip
- HU** Többlumenű, centrális vénás katéterezési szett, Blue FlexTip katéterrel
- IS** Fjöllumena holæðaleggsett með Blue FlexTip æðalegg
- IT** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- LT** Daugiakanalis centrinės venos kateterizavimo mini rinkinys su „Blue FlexTip“ kateteriu
- NO** Multilumen sentralt venekateriseringssett med Blue FlexTip kateter
- PL** Zestaw do cewnikowania żyły centralnej z cewnikiem wielokanałowym z końcówką Blue FlexTip
- PT** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip
- RO** Set de cateterizare venoasă centrală multi-lumen cu cateter Blue FlexTip
- RU** Набор с многосветным центральным венозным катетером с кончиком Blue FlexTip
- SR** Višelumenski centralni venski kateterizacijski set sa Blue FlexTip kateterom
- SK** Vialúmenová centrálna venózna katetrizačná súprava s katetrom Blue FlexTip
- SL** Večlumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip
- ES** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- SV** Set för flerkanalig central venkatetrering med Blue FlexTip kateter
- TR** Multi Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- UK** Багатоприсвітний центральний венозний катетер Набор: катетер Blue с кінчиком Blue FlexTip Catheter



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)
Distal (16 Ga)	0.47	3530
Medial (14 Ga)	0.89	7555
Proximal (16 Ga)	0.47	3290

\* Priming volumes are approximate and are done without the injection cap.  
† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

57-15853-115R (11/14)

