

Dimensions:
Width: 8-3/4" (22.2 cm)
Length: 11-1/8" (28.3 cm)

REF EU-15853-CVET

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

ARROW®

Contents:

- 1: Multi-Lumen Indwelling Catheter:
8.5 Fr x 8" (20 cm) Radiopaque
Polyurethane with Blue FlexTip®,
Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032"
(0.81 mm) dia. x 23-5/8" (60 cm)
(Straight Soft Tip on One End — "J" Tip on
Other) with Arrow Advance™
with ECG Mark
- 1: Introducer Needle: Echogenic 18 Ga
x 2-1/2" (6.35 cm) XTW
- 1: Arrow® Raulerson Spring-Wire
Introduction Syringe: 5 ml
- 1: Tissue Dilator 9 Fr (3.0 mm) x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Scalpel: #11
- 1: ECG Extension Cable
- 1: Pressure Transduction Probe
- 2: Injection Site Cap

All components are CE 0086 unless otherwise noted.

Fluid path components are non-pyrogenic.




Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Teleflex®
www.teleflex.com

- be** Многопросветный набор для центральной венозной катетеризации с катетером Blue FlexTip
- hr** Višelumenski centralni venski kateter Set s Blue FlexTip kateterom
- cs** Vícelumenový centrální žilní katetizační set s katetrem Blue FlexTip
- da** Multilumen centralvenøs kateteriseringsset med Blue FlexTip kateter
- nl** Meerlumen centraal veneuze katheterisatieset met Blue FlexTip katheter
- et** Mitme valendikuga tsentraalse veenikateeteri komplekt koos Blue FlexTip kateetriga
- fi** Monilumeninen keskuslaskimokatetri sarja jossa Blue FlexTip katetri
- fr** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- de** Besteck für viellumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- el** Σετ Καθητηριασμού Κεντρικής φλέβας Πολλαπλούς Αυλούς με Καθητήρα Blue FlexTip
- hu** Többlumenű, centrális vénás katéterezési szett, Blue FlexTip katéterrel
- is** Fjöl lúmena holæðaleggset með Blue FlexTip æðalegg
- it** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- lt** Daugiakanalis centrinės venos kateterizavimo mini rinkinys su „Blue FlexTip“ kateteriu
- no** Multilumen sentralt venekateteriseringssett med Blue FlexTip kateter
- pl** Zestaw do cewnikowania żyły centralnej z cewnikiem wielokanałowym z końcówką Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip
- ro** Set de cateterizare venoasă centrală multi-lumen cu cateter Blue FlexTip
- ru** Набор с многопросветным центральным венозным катетером с кончиком Blue FlexTip
- sr** Višelumenski centralni venski kateterizacioni Set sa Blue FlexTip kateterom
- sk** Viačlúmenová centrálna venózna katetizačná súprava s katétrom Blue FlexTip
- sl** Večlumenski centralni venski komplet za katetizacijo s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- sv** Set för flerkanalig central venkateterisering med Blue FlexTip kateter
- tr** Multi Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- uk** Багатоприсвітний центральний венозний катетер Набор: катетер Blue с кінчиком Blue FlexTip Catheter



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga)	0.47	3530
Medial (14 Ga)	0.89	7555
Proximal (16 Ga)	0.47	3290

* Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.
† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



52-15853-1188 (11/16)