

Dimensions:
Width: 8-3/4" (22.2 cm)
Length: 11-1/8" (28.3 cm)

REF EU-15854-CVET
Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter
ARROW®

Contents:

- 1: Quad-Lumen Indwelling Catheter: 8.5 Fr x 8" (20 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps, Injection Site Caps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™ with ECG Mark
- 1: Introducer Needle: Echogenic 18 Ga x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 ml
- 1: Tissue Dilator 9 Fr (3.0 mm) x 4" (10.2 cm)
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter
- 1: ECG Extension Cable
- 1: Scalpel: #11

All components are CE 0086 unless otherwise noted.


Contains no medication.
Fluid path components are non-pyrogenic.



Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

- BE** Чатырохпросветны набор для цэнтральнай венознай катэтарызацыі з катэтрам Blue FlexTip
- HR** Quad-Lumen centralni venski kateter Komplet s Blue FlexTip kateter
- CS** Čtyřlumenový centrální žilní katetrizační set s katetrem Blue FlexTip
- DA** Firdobbelt centralvenøs kateteriserings-sæt med Blue FlexTip kateter
- NL** Centraal veneuze quad-lumen katheterisatieset met Blue FlexTip katheter
- ET** Neljavalendikulise tsentraalveeni kateteriseerimiskomplekt Blue FlexTip kateetriga
- F** Quad-Lumen keskuslaskimokatetri sarja Blue FlexTip katetrilla
- FR** Jeu de cathétérisme veineux central à quatre lumières avec cathéter Blue FlexTip
- DE** Besteck für vierlumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- EL** Σετ καθετηριασμού κεντρικής φλέβας τετραπλάσιο αυλού με καθετήρα Blue FlexTip
- HU** Quad-Lumen centrális vénás katéterezési szett, Blue FlexTip katéterrel
- IS** Fjögur lúmena holæðaleggsett með Blue FlexTip æðalegg
- IT** Set per cateterismo della vena centrale a quattro lumi con catetere Blue FlexTip
- LT** Keturių angų centrinės venos kateterizavimo rinkinys su „Blue FlexTip“ kateteriu
- NO** Sentralt venekateteriseringssett med fire lumen og Blue FlexTip kateter
- PL** Zestaw do cewnikowania żył centralnych z 4-kanalowym cewnikiem z końcówką Blue FlexTip
- PT** Conjunto de Cateterização Venosa Central com Lúmen Quádruplo e Cateter Blue FlexTip
- RO** Set de cateterizare venoasă centrală Quad-Lumen cu cateter Blue FlexTip
- RU** Четырех-просветный комплект для установки центрального венозного катетера с катетером Blue FlexTip
- SR** Quad-Lumen komplet za centralnu vensku kateterizaciju sa Blue FlexTip kateterom
- SK** Štvorlúmenová centrálna venózna katetrizačná súprava s katétrom Blue FlexTip
- SL** Štiri-lumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip
- ES** Conjunto de cateterización venosa central de cuatro luces con catéter Blue FlexTip
- SV** Set för fyrkanalyg central venkateterisering med Blue FlexTip kateter
- TR** Dörtlü-Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- UK** Чотирипросвітний комплект для встановлення центрального венозного катетера з катетером Blue FlexTip



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)
Distal (16 Ga)	0.43	2700
Medial 1 (14 Ga)	0.54	5100
Medial 2 (18 Ga)	0.36	1400
Proximal (18 Ga)	0.41	1600

* Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.
† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

LOT



(01) 0 08 01902 07627 4

SV15854-17A (4/14)



REF EU-15854-CVET 8.5 Fr 4 Lumen 20 cm catheter length .032 inch dia. spring-wire guide
Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open