

Dimensions:
Width: 6-3/4" (17.2 cm)
Length: 11-5/16" (28.7 cm)

REF EU-15955-N
Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter
ARROW®

Contents:

- 1: Multi-Lumen Indwelling Catheter: 9.5 Fr x 8" (20 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™ with ECG Mark
- 1: Introducer Needle: 18 Ga x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 ml Luer-Slip
- 1: Tissue Dilator 10 Fr (3.5 mm) x 4" (10.2 cm)
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter

All components are CE 0086 unless otherwise noted.


Contains no medication.
Fluid path components are non-pyrogenic.



Arrow International, Inc.
2400 Bemville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

- BE** Многопросветный набор для центральной венозной катетеризации з катэтрам Blue FlexTip
- HR** Višelumenski centralni venski kateter Set s Blue FlexTip kateterom
- CS** Vícelumenový centrální žilní katetrizační set s katetrem Blue FlexTip
- DA** Multilumen centralvenøs kateteriseringsset med Blue FlexTip kateter
- NL** Multi-lumen centraal veneuze katheterisatieset met Blue FlexTip katheter
- ET** Mitme valendikuga tsentraalse veenikateeteri komplekt koos Blue FlexTip kateetriga
- F** Monilumeninen keskuskaskimokatetri sarja jossa sininen FlexTip katetri
- FR** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- DE** Besteck für viellumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- EL** Σετ Καθετηριασμού Κεντρικής φλέβας Πολλαπλούς Αυλούς με Καθετήρα Blue FlexTip
- HU** Többlumenű, centrális vénás katéterezési szett, Blue FlexTip katéterrel
- IS** Fjöllumena holæðaleggsett með Blue FlexTip æðalegg
- IT** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- LT** Daugiakanaliscentrinės venos kateterizavimo mini rinkinys su „Blue FlexTip“ kateteriu
- NO** Multilumen sentralt venekateriseringssett med Blue FlexTip kateter
- PL** Zestaw do cewnikowania żyły centralnej z cewnikiem wielokanałowym z końcówką Blue FlexTip
- PT** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip
- RO** Set de cateterizare venoasă centrală multi-lumen cu cateter Blue FlexTip
- RU** Набор с многопросветным центральным венозным катетером с кончиком Blue FlexTip
- SR** Višelumenski centralni venski kateterizacijski Set sa Blue FlexTip kateterom
- SK** Vialúmenová centrálna venózna katetrizačná súprava s katétrom Blue FlexTip
- SL** Večlumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip
- ES** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- SV** Set för flerkanalig central venkateterisering med Blue FlexTip kateter
- TR** Multi Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateteri ile
- UK** Багатопроевітний центральний венозний катетер Набір: катетер Blue с кінчиком Blue FlexTip Catheter



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)
Distal (16 Ga)	0.48	3605
Medial 1 (14 Ga)	0.57	5390
Medial 2 (18 Ga)	0.38	1250
Medial 3 (18 Ga)	0.38	1160
Proximal (18 Ga)	0.41	1455

*Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.20 ml.
†Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

SZ-15955-119A (4/14)



REF EU-15955-N 9.5 Fr 5 Lumen 20 cm catheter length .032 inch dia. spring-wire guide
Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open